


DESS














> DENTAL SMART SOLUTIONS




Catalogue 2022/2

PRODUCT LIST



 NOBELACTIVE® & NOBELREPLACE® CC ACTIVE HEX 04	 ZIMMER SCREW-VENT® INTERNAL HEX USA 54	 BIOHORIZONS® EXTERNAL EXTERNAL HEX BH 98	 BEGO SEMADOS® INTERNAL SEM 128
 NOBEL REPLACE SELECT™ TRI-LOBE 10	 ASTRA TECH IMPLANT SYSTEM™ EV CONIC EVO 60	 MIS® SEVEN INTERNAL HEX MI 102	 BIOTECH KONTACT® CONIC BIO 130
 NOBEL BRÅNEMARK® EXTERNAL HEX UNIVERSAL 16	 ASTRA TECH OSSEOSPEED™ INTERNAL HEX CONIC 66	 MIS® C1 CONIC MC 108	 BTI® INTERNAL BT 132
 3i OSSEOTITE® EXTERNAL HEX USA 22	 DENTSPLY FRIADENT XIVE® INTERNAL HEX FD 72	 MIS® V3 CONIC MV 110	 PHIBO® TSH® EXTERNAL HEX PH 134
 3i CERTAIN® INTERNAL HEX "CLICK" 28	 DENTSPLY ANKYLOS® C/X INTERNAL ANK 78	 MEGAGEN AnyRidge® CONIC ANYR 112	 MULTI-UNIT CONNECTION 136
 STRAUMANN® SOFT TISSUE LEVEL & synOcta® OCTAGON 34	 OSSTEM® TS CONIC OSS 82	 ANTHOGYR AXIOM® CONIC AXIO 116	 DESSLoc® 140
 STRAUMANN® BONE LEVEL CONICAL BL 40	 CAMLOG® INTERNAL CAM 86	 MEDENTIS ICX® CONIC IC 118	 HOLDERS & BLOOCKS 144
 STRAUMANN® BLX CONICAL BLX 46	 CONELOG® INTERNAL CON 90	 GLOBALD IN-KONE® CONIC INK 122	 TOOLS 146
 NEODENT® GRAND MORSE NEO GM 50	 BIOHORIZONS® INTERNAL INTERNAL HEX BH 94	 DENTIUM SUPERLINE™ & IMPLANTIUM® CONIC DENT 124	 CLINICAL SCREWS RECOMMENDED TORQUE 150



ACTIVE HEX


compatible with
NobelActive® / Replace® CC

- non-engaging
- engaging



HEALING ABUTMENT	3.0	NP	RP	WP
	3.0	3.5	4.3/5.0	5.5
GH 3,0mm ^A	20.068/3	20.039/3	20.040/3	20.069/3
GH 3,0mm ^B		20.047/3	20.048/3	20.070/3
GH 4,0mm ^A	20.068/4	20.039/4	20.040/4	20.069/4
GH 4,0mm ^B		20.047/4	20.048/4	20.070/4
GH 5,0mm ^A	20.068/5	20.039/5	20.040/5	20.069/5
GH 5,0mm ^B		20.047/5	20.048/5	20.070/5

^A for single crown restoration ^B for bridges and structures restoration



SCAN ABUTMENTS					
		50.068	50.041	50.042	50.069
		52.068	52.041	52.042	52.069
		53.100-P10	53.100-P10	53.100-P10	53.100-P10
		53.200-P10	53.200-P10	53.200-P10	53.200-P10



TRANSFER				
	17.068	17.044	17.045	17.069

Transfers include two screws (for open and closed tray).



Platform Colour Coded

DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59



Platform Colour Coded

ANALOGUE

3.0	NP	RP	WP
3.0	3.5	4.3/5.0	5.5
14.068/D	14.041/D	14.042/D	14.069/D
10 Pack.			
	14.041-P10	14.042-P10	14.069-P10



TEMPORARY ABUTMENT

	24.041	24.042
	25.041	25.042



DESSLoc®


GH 1,0mm	88.039/1	88.040/1
GH 2,0mm	88.039/2	88.040/2
GH 3,0mm	88.039/3	88.040/3
GH 4,0mm	88.039/4	88.040/4



MULTI-UNIT ABUTMENTS			3.0	NP	RP	WP
			3.0	3.5	4.3/5.0	5.5
GH 1,5mm	0°			40.007/1	40.008/1	
GH 2,5mm	0°			40.007/2	40.008/2	
GH 3,5mm	0°			40.007/3	40.008/3	
GH 4,5mm	0°				40.008/4	
GH 2,5mm	17°			41.707/2	41.708/2	
GH 3,5mm	17°			41.707/3	41.708/3	
GH 3,5mm	30°			43.007/3	43.008/3	
GH 4,5mm	30°			43.007/4	43.008/4	
			4.8 MULTI-UNIT PLATFORM			

TI-BASE

GH 1,2mm			15.041	15.042	
GH 1,5mm		15.068			15.069
GH 2,0mm			15.041/2	15.042/2	
GH 3,0mm			15.341	15.342	
GH 1,2mm			16.041	16.042	
GH 1,5mm		16.068			16.069
GH 2,0mm			16.041/2	16.042/2	
GH 3,0mm			16.341	16.342	
TI-BASE CASTABLES					
5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5




AURUMBase®

		3.0	NP	RP	WP
		3.0	3.5	4.3/5.0	5.5
			36.041	36.042	
			37.041	37.042	
PLASTIC SLEEVES FOR AURUMBase®					
 5 Pack.	0°		33.100-P5	33.100-P5	
 5 Pack.	10°		33.101-P5	33.101-P5	
 5 Pack.	20°		33.102-P5	33.102-P5	

ELLIPTIBase®

A corner		56.068
B flat		57.068



C-Base®

GH 0,8mm		58.041	58.042
GH 3,0mm		58.041/3	58.042/3
GH 0,8mm		59.041	59.042
GH 3,0mm		59.041/3	59.042/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.

Torx® Screw
AURUMBase® & ELLIPTIBase®

(replacement screw)	19.468	19.441	19.442
Recommended torque	35 Ncm	35 Ncm	35 Ncm



Torx® BALL SCREWDRIVER

	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06





ABUTMENTS	3.0	NP	RP	WP
	3.0	3.5	4.3/5.0	5.5
	0°	13.041	13.042	
	15°	22.041	22.042	
	25°		23.042	



SCREW UG	19.068	19.027	19.030	19.030
	10 Pack.	19.027-P10	19.030-P10	19.030-P10
DLC	19.268	19.227	19.230	19.230
TIN	19.668	19.627	19.630	19.630
Recommended torque	15 Ncm	35 Ncm	35 Ncm	35 Ncm



PRE-MILLED BLANK					
Ø 10mm	Ti	61.068	61.041	61.042	61.069
Ø 14mm	Ti		62.041	62.042	62.069
Ø 14mm	CoCr	71.068	71.041	71.042	71.069



CoCr BASE				
	○	30.041	30.042	
	⬡	31.041	31.042	



CASTABLE				
	○	28.041	28.042	
		11.041-P10	11.042-P10	
	⬡	29.041	29.042	
		12.041-P10	12.042-P10	


“Staying true to our history and values empowers our brand”

Since our beginnings, quality has been quintessential to our philosophy and is reflected in every product that we manufacture. Our choice of the best Swiss and Japanese milling machines allow us to create the most precise parts. Our quality standards entail measuring to within a few microns of precision.





HEALING ABUTMENT	NP	RP	WP	6.0
	NP/3.5	RP/4.3	WP/5.0	6.0
GH 3,0mm	20.004/3	20.005/3	20.006/3	
GH 4,0mm	20.004/4	20.005/4	20.006/4	
GH 5,0mm	20.004/5	20.005/5	20.006/5	



SCAN ABUTMENTS						
			50.004	50.005	50.006	50.097
			52.004	52.005	52.006	52.097
	10 Pack.		53.100-P10	53.100-P10	53.100-P10	53.100-P10
	10 Pack.		53.200-P10	53.200-P10	53.200-P10	53.200-P10




TRANSFER				
	17.004	17.005	17.006	17.097
<i>Transfers include two screws (for open and closed tray).</i>				



	Platform Colour Coded	Retention screw Hex 1,27 mm*	14.004/D	14.005/D	14.006/D	14.097/D
*Other retention screws available, see page 59						



ANALOGUE	NP	RP	WP	6.0
	NP/3.5	RP/4.3	WP/5.0	6.0
	14.004	14.005	14.006	14.097
10 Pack.	14.004-P10	14.005-P10		



TEMPORARY ABUTMENT				
	24.004	24.005	24.006	
	25.004	25.005	25.006	



DESSLoc®		
GH 1,0mm	88.004/1	88.005/1
GH 2,0mm	88.004/2	88.005/2
GH 3,0mm	88.004/3	88.005/3
GH 4,0mm	88.004/4	88.005/4



MULTI-UNIT
ABUTMENTS

			NP		RP		WP		6.0
			NP/3.5		RP/4.3		WP/5.0		6.0
GH 1,0mm	0°		40.004/1		40.005/1		40.006/1		
GH 2,0mm	0°		40.004/2		40.005/2		40.006/2		
GH 3,0mm	0°		40.004/3		40.005/3		40.006/3		
GH 4,0mm	0°				40.005/4				
GH 5,0mm	0°				40.005/5				
GH 2,0mm	17°		41.704/2		41.705/2				
GH 3,0mm	17°		41.704/3		41.705/3				
GH 4,0mm	17°				41.705/4				
GH 4,0mm	30°				43.005/4				
GH 5,0mm	30°				43.005/5				
GH 5,0mm	30°				43.005/5R				



4.8
MULTI-UNIT PLATFORM

TI-BASE

GH 0,3mm		15.004	15.005	15.006	15.097
GH 1,5mm		15.204	15.205	15.206	
GH 3,0mm		15.304	15.305	15.306	
GH 0,3mm		16.004	16.005	16.006	16.097
GH 1,5mm		16.204	16.205	16.206	
GH 3,0mm		16.304	16.305	16.306	

TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5



AURUMBase®

		NP		RP		WP		6.0
		NP/3.5		RP/4.3		WP/5.0		6.0
		36.004		36.005		36.006		
		37.004		37.005		37.006		

PLASTIC SLEEVES FOR AURUMBase®

	5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5

C-Base®

GH 0,4mm		58.004		58.006	58.097
GH 0,5mm			58.005		
GH 1,5mm		58.004/2	58.005/2	58.006/2	
GH 3,0mm		58.004/3	58.005/3	58.006/3	
GH 0,4mm		59.004		59.006	59.097
GH 0,5mm			59.005		
GH 1,5mm		59.004/2	59.005/2	59.006/2	
GH 3,0mm		59.004/3	59.005/3	59.006/3	

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.

Torx® Screw AURUMBase®

(replacement screw)	19.404	19.405	19.405
Recommended torque	35 Ncm	35 Ncm	35 Ncm

Torx® BALL SCREWDRIVER


20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS		NP	RP	WP	6.0
		NP/3.5	RP/4.3	WP/5.0	6.0
	0°	13.004	13.005	13.006	
	15°	22.004	22.005	22.006	
	25°		23.005		



SCREW UG					
		19.004	19.005	19.005	19.005
	10 Pack.	19.004-P10	19.005-P10	19.005-P10	19.005-P10
	DLC	19.204	19.205	19.205	19.205
	TIN	19.604	19.605	19.605	19.605
Recommended torque		35 Ncm	35 Ncm	35 Ncm	35 Ncm



PRE-MILLED BLANK					
Ø 10mm	Ti	61.004	61.005	61.006	61.097
Ø 14mm	Ti		62.005	62.006	62.097
Ø 14mm	CoCr	71.004	71.005	71.006	71.097



CoCr BASE					
		30.004	30.005	30.006	
		31.004	31.005	31.006	



CASTABLE					
		28.004	28.005	28.006	
		11.004-P10	11.005-P10		
		29.004	29.005	29.006	
		12.004-P10	12.005-P10		

“Great things are done by a series of small things brought together”

A business without quality will not last. Terrats Medical is in its third generation and we want to honour our heritage with the highest requirements for quality in all we do.





EXTERNAL HEX UNIVERSAL


compatible with
Nobel Brånemark®

- non-engaging
- ⬡ engaging



HEALING ABUTMENT

	NP	RP	WP
	NP/3.5	RP/4.1	WP/5.1
GH 3,0mm	20.001/3	20.002/3	20.003/3
GH 4,0mm	20.001/4	20.002/4	20.003/4
GH 5,0mm	20.001/5	20.002/5	20.003/5



SCAN ABUTMENTS

		50.001	50.002	50.003
		52.001	52.002	52.003
	10 Pack. ○	53.100-P10	53.100-P10	53.100-P10
	10 Pack. ⬡	53.200-P10	53.200-P10	53.200-P10



TRANSFER

	17.041	17.042	17.043
--	--------	--------	--------

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE

	NP	RP	WP
	NP/3.5	RP/4.1	WP/5.1
Platform Colour Coded	14.001/D	14.002/D	14.003/D

*Other retention screws available, see page 59



ANALOGUE

Platform Colour Coded	14.001	14.002	14.003
10 Pack.		14.002-P10	




TEMPORARY ABUTMENT

○	24.001	24.002	24.003
⬡	25.001	25.002	25.003

DESSLoc®

GH 1,0mm	88.001/1	88.002/1
GH 2,0mm	88.001/2	88.002/2
GH 3,0mm	88.001/3	88.002/3
GH 4,0mm		88.002/4
GH 5,0mm		88.002/5





MULTI-UNIT ABUTMENTS				NP	RP	WP
				NP/3.5	RP/4.1	WP/5.1
GH 1,0mm	0°	○		40.001/1	40.002/1	40.003/1
GH 2,0mm	0°	○		40.001/2	40.002/2	40.003/2
GH 3,0mm	0°	○		40.001/3	40.002/3	40.003/3
GH 4,0mm	0°	○			40.002/4	
GH 5,0mm	0°	○			40.002/5	
GH 2,0mm	17°	⬡		41.701/2	41.702/2	
GH 3,0mm	17°	⬡		41.701/3	41.702/3	
GH 4,0mm	17°	⬡			41.702/4	
GH 4,0mm	30°	⬡			43.002/4	
GH 5,0mm	30°	⬡			43.002/5	

4.8

6.5


MULTI-UNIT PLATFORM



TI-BASE						
GH 0,3mm		○		15.001	15.002	15.003
GH 1,5mm		○		15.201	15.202	15.203
GH 3,0mm		○		15.301	15.302	15.303
GH 0,3mm		⬡		16.001	16.002	16.003
GH 1,5mm		⬡		16.201	16.202	16.203
GH 3,0mm		⬡		16.301	16.302	16.303

TI-BASE CASTABLES

5 Pack.	○	32.001-P5	32.001-P5	32.001-P5
5 Pack.	⬡	33.001-P5	33.001-P5	33.001-P5



AURUMBase®			NP	RP	WP
			NP/3.5	RP/4.1	WP/5.1
GH 0,3mm	○		36.001	36.002	
GH 1,5mm	○		36.001/2	36.002/2	
GH 2,5mm	○		36.001/3	36.002/3	
GH 0,3mm	⬡		37.001	37.002	
GH 1,5mm	⬡		37.001/2	37.002/2	
GH 2,5mm	⬡		37.001/3	37.002/3	

PLASTIC SLEEVES FOR AURUMBase®


5 Pack.	0°	33.100-P5	33.100-P5
5 Pack.	10°	33.101-P5	33.101-P5
5 Pack.	20°	33.102-P5	33.102-P5



C-Base®			
GH 1,5mm	○	58.001/2	58.002/2
GH 3,0mm	○	58.001/3	58.002/3
GH 0,6mm	⬡	59.001	59.002
GH 1,5mm	⬡	59.001/2	59.002/2
GH 3,0mm	⬡	59.001/3	59.002/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



Torx® Screw AURUMBase®		
(replacement screw)	19.401	19.402
Recommended torque	35 Ncm	35 Ncm

Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06





ABUTMENTS

		NP	RP	WP
		NP/3.5	RP/4.1	WP/5.1
	0°	13.001	13.002	13.003
	15°	22.001	22.002	22.003
	25°		23.002	



SCREW UG

		19.027	19.028	19.029
10 Pack.		19.027-P10	19.028-P10	19.029-P10
DLC		19.227	19.228	19.229
TIN		19.627	19.628	19.629
Recommended torque		35 Ncm	35 Ncm	35 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.001	61.002	61.003
Ø 14mm	CoCr	71.001	71.002	71.003



CoCr BASE

	○	30.001	30.002	30.003
	⬡	31.001	31.002	31.003



FOR WELDING TECHNIQUES

		○	30.102	
--	--	---	--------	--




CASTABLE

10 Pack.	○	28.001	28.002	28.003
			11.002-P10	
10 Pack.	⬡	29.001	29.002	29.003
			12.002-P10	

“The difference between something good and something exceptional is in the smallest detail”


At DESS® we believe that no detail is too small to escape close attention. That’s why our exhaustive quality controls in every single part of the manufacturing process ensure that the quality of our products is guaranteed. Our confidence is such that we offer a lifetime warranty on all our products.





HEALING ABUTMENT

	NP	RP	WP
	3.4	4.1	5.0
GH 3,0mm	20.011/3	20.012/3	20.013/3
GH 4,0mm	20.011/4	20.012/4	20.013/4
GH 5,0mm	20.011/5	20.012/5	20.013/5



SCAN ABUTMENTS


		50.011	50.012	50.013
		52.011	52.012	52.013
	10 Pack. 	53.100-P10	53.100-P10	53.100-P10
	10 Pack. 	53.200-P10	53.200-P10	53.200-P10



TRANSFER


17.011	17.012	17.013
--------	--------	--------

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*
*Other retention screws available, see page 59



ANALOGUE


Platform Colour Coded

10 Pack.



TEMPORARY ABUTMENT

	24.011	24.012	24.013
	25.011	25.012	25.013



DESSLoc®

GH 1,0mm	88.011/1	88.002/1	88.013/1
GH 2,0mm	88.011/2	88.002/2	88.013/2
GH 3,0mm	88.011/3	88.002/3	88.013/3
GH 4,0mm	88.011/4	88.002/4	88.013/4
GH 5,0mm	88.011/5	88.002/5	



MULTI-UNIT
ABUTMENTS

			NP	RP	WP
			3.4	4.1	5.0
GH 1,0mm	0°		40.009/1	40.017/1	
GH 2,0mm	0°		40.009/2	40.017/2	
GH 3,0mm	0°		40.009/3	40.017/3	
GH 4,0mm	0°			40.017/4	
GH 5,0mm	0°			40.017/5	
GH 2,0mm	17°			41.712/2	
GH 3,0mm	17°			41.712/3	
GH 4,0mm	17°			41.712/4	
GH 3,0mm	30°			43.012/3	
GH 4,0mm	30°			43.012/4	
GH 5,0mm	30°			43.012/5	



4.8
MULTI-UNIT PLATFORM

TI-BASE

GH 0,3mm		15.011	15.012	15.013
GH 1,5mm		15.211	15.212	15.213
GH 3,0mm		15.311	15.312	15.313
GH 0,3mm		16.011	16.012	16.013
GH 1,5mm		16.211	16.212	16.213
GH 3,0mm		16.311	16.312	16.313

TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5



AURUMBase®

		NP	RP	WP
		3.4	4.1	5.0
GH 0,3mm		36.011	36.012	36.013
GH 1,5mm		36.011/2	36.012/2	
GH 2,5mm		36.011/3	36.012/3	
GH 0,3mm		37.011	37.012	37.013
GH 1,5mm		37.011/2	37.012/2	
GH 2,5mm		37.011/3	37.012/3	

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5

C-Base®

GH 1,5mm		58.011/2	58.012/2	58.013/2
GH 3,0mm		58.011/3	58.012/3	58.013/3
GH 0,6mm		59.011	59.012	59.013
GH 1,5mm		59.011/2	59.012/2	59.013/2
GH 3,0mm		59.011/3	59.012/3	59.013/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.

Torx® Screw AURUMBase®

(replacement screw)	19.402	19.402	19.402
Recommended torque	20 Ncm	20 Ncm	20 Ncm

Torx® BALL SCREWDRIVER

20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS

		NP	RP	WP
		3.4	4.1	5.0
	0°	13.011	13.012	13.013
	15°	22.011	22.012	22.013
	25°		23.012	



SCREW HEX 1,20mm

		19.034	19.034	19.034
10 Pack.		19.034-P10	19.034-P10	19.034-P10
TIN		19.634	19.634	19.634
Recommended torque		20 Ncm	20 Ncm	20 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.011	61.012	61.013
Ø 14mm	CoCr	71.011	71.012	71.013



CoCr BASE

	○	30.011	30.012	30.013
	⬡	31.011	31.012	31.013



CASTABLE

10 Pack.	○	28.011	28.012	28.013
			11.002-P10	
10 Pack.	⬡	29.011	29.012	29.013
			12.002-P10	

“There is no good idea that cannot be improved”

We are proud to say that innovation and quality are at the center of our core values. We look at different solutions and re-engineer them to provide better and more reliable alternatives. Products like our intra-oral scan abutment without a screw hole, which improves substantially the image precision, our angled Multi-unit abutment with a fully shaped connection and vertical access for more stability and less screw loosening or our angulated screw channel solution AURUMBase® that allows you 25° angle with a 360° rotation, are clear examples of our way of doing things.





EXTERNAL HEX “CLICK”

compatible with
3i Certain®

-  non-engaging
-  engaging



HEALING ABUTMENT

	NP	RP	WP
	3.4	4.1	5.0
GH 3,0mm	20.014/3	20.015/3	20.016/3
GH 4,0mm	20.014/4	20.015/4	20.016/4
GH 5,0mm	20.014/5	20.015/5	20.016/5



SCAN ABUTMENTS



		50.014	50.015	50.016
		52.014	52.015	52.016
	10 Pack. 	53.100-P10	53.100-P10	53.100-P10
	10 Pack. 	53.200-P10	53.200-P10	53.200-P10



TRANSFER

	17.014	17.015	17.016
--	--------	--------	--------

Transfers include two screws (for open and closed tray).



Platform Colour Coded

DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59



Platform Colour Coded

ANALOGUE



TEMPORARY ABUTMENT



	24.014	24.015	24.016
	25.014	25.015	25.016




DESSLoc®


GH 1,0mm	88.014/1	88.015/1	88.016/1
GH 2,0mm	88.014/2	88.015/2	88.016/2
GH 3,0mm	88.014/3	88.015/3	88.016/3
GH 4,0mm	88.014/4	88.015/4	88.016/4
GH 5,0mm	88.014/5	88.015/5	



			NP	RP	WP
			3.4	4.1	5.0
TI-BASE					
	GH 0,3mm	○	15.014	15.015	15.016
	GH 1,5mm	○	15.214	15.215	15.216
	GH 3,0mm	○	15.314	15.315	15.316
	GH 0,3mm	⬡	16.014	16.015	16.016
	GH 1,5mm	⬡	16.214	16.215	16.216
	GH 3,0mm	⬡	16.314	16.315	16.316
TI-BASE CASTABLES					
	5 Pack.	○	32.001-P5	32.001-P5	32.001-P5
	5 Pack.	⬡	33.001-P5	33.001-P5	33.001-P5



AURUMBase®					
		○	36.014	36.015	36.016
		⬡	37.014	37.015	37.016
PLASTIC SLEEVES FOR AURUMBase®					
	5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5



			NP	RP	WP
			3.4	4.1	5.0
C-Base®					
	GH 0,4mm	○	58.014	58.015	
	GH 0,5mm	○			58.016
	GH 1,5mm	○	58.014/2	58.015/2	58.016/2
	GH 3,0mm	○	58.014/3	58.015/3	58.016/3
	GH 0,4mm	⬡	59.014	59.015	
	GH 0,5mm	⬡			59.016
	GH 1,5mm	⬡	59.014/2	59.015/2	59.016/2
	GH 3,0mm	⬡	59.014/3	59.015/3	59.016/3
Includes screw, same as standard screw.			Libraries available for EXOCAD & 3Shape.		



Torx® Screw AURUMBase®			
(replacement screw)	19.438	19.438	19.438
Recommended torque	20 Ncm	20 Ncm	20 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS				
	0°	13.014	13.015	13.016
	15°	22.014	22.015	22.016
	25°		23.015	



SCREW HEX 1,20mm

	NP	RP	WP
	3.4	4.1	5.0
	19.011	19.011	19.011
10 Pack.	19.011-P10	19.011-P10	19.011-P10
TIN	19.611	19.611	19.611
Recommended torque	20 Ncm	20 Ncm	20 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.014	61.015	61.016
Ø 14mm	Ti		62.015	62.016
Ø 14mm	CoCr	71.014	71.015	71.016



CoCr BASE

	30.014	30.015	30.016
	31.014	31.015	31.016



CASTABLE

	28.014	28.015	28.016
10 Pack.		11.015-P10	
	29.014	29.015	29.016
10 Pack.		12.015-P10	


Your design, our precision

Mill your design with the precision of the DESS® connections. Our range of blank premilled is available for the most popular implant connections in the world.



- 32 compatibilities
- Titanium Grade 5 ELI Ø 10mm 14mm
- Micromelt CoCr Ø 14mm

- Total millable lenght 20mm
- Minimal gingival transition (0,3mm)



OCTAGON

compatible with
Straumann® Tissue level & synOcta®

- non-engaging
- engaging




HEALING
ABUTMENT

	NNC	RN	WN
	3.5	4.1/4.8	4.8/6.5
GH 2,0mm		20.009/2	20.045/2
GH 3,0mm		20.009/3	20.045/3
GH 4,0mm		20.009/4	20.045/4
On abutment		20.010	20.049



CONVERTER
ABUTMENT

		35.001	
		34.001	34.002



SCAN
ABUTMENTS

		50.070	50.009	50.045
On abutment			50.010	50.048
		52.070	52.009	52.045
On abutment			52.010	52.048
	10/5 Pack.	53.100-P10	53.110-P5	53.148-P5
	10 5 Pack.	53.200-P10	53.210-P5	53.248-P5



TRANSFER

	NNC	RN	WN
	3.5	4.1/4.8	4.8/6.5
	17.070	17.009	17.053
On abutment		17.010	17.048

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*	14.070/D	14.009/D	14.045/D
On abutment		14.010/D	14.046/D

*Other retention screws available, see page 59



ANALOGUE




Platform Colour Coded

	14.070	14.009	14.045
	10 Pack.	14.009-P10	
On abutment		14.010	14.046
	10 Pack.	14.010-P10	



TEMPORARY
ABUTMENT


		24.026	24.045
		25.009	25.045
On abutment		24.010	24.047
		25.026	25.047



DESSLoc®

	NNC	RN	WN
	3.5	4.1/4.8	4.8/6.5
GH 1,0mm		88.009/1	
GH 2,0mm		88.009/2	
GH 3,0mm		88.009/3	
GH 4,0mm		88.009/4	

TI-BASE




GH 0,3mm	○		15.045
GH 0,4mm	○		15.046
GH 0,5mm	○	15.070	
GH 0,3mm	⬡		16.045
GH 0,4mm	⬡		16.009
GH 0,5mm	⬡	16.070	
On abutment GH 0,3mm	○	15.010	15.048
GH 0,3mm	⬡	16.010	16.048



TI-BASE CASTABLES

5 Pack.	○	32.001-P5	32.046-P5	32.045-P5
5 Pack.	⬡	33.001-P5	33.009-P5	33.045-P5
On abutment 5 Pack.	○		32.010-P5	32.045-P5
5 Pack.	⬡		33.010-P5	33.045-P5




AURUMBase®

	○	36.046	36.047
	⬡	37.009	37.045

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.200-P5	33.300-P5
5 Pack.	10°	33.201-P5	33.301-P5
5 Pack.	20°	33.202-P5	33.302-P5




C-Base®

	NNC	RN	WN
	3.5	4.1/4.8	4.8/6.5
GH 1,0mm	○	58.009	
GH 1,2mm	○		58.045
GH 1,0mm	⬡	59.009	
GH 1,2mm	⬡		59.045

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



A

B

ELLIPTIBase®

A corner	⬡	56.070
B flat	⬡	57.070



Torx® Screw
AURUMBase® & ELLIPTIBase®

(replacement screw)	19.441	19.446	19.446
Recommended torque	35 Ncm	35 Ncm	36 Ncm



Torx® BALL SCREWDRIVER

20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS

		NNC	RN	WN
		3.5	4.1/4.8	4.8/6.5
	0°		13.009*	
	15°		22.009*	
	25°		23.009*	
		*inc.captive screw		



SCREW TORX®

		19.070	19.026	19.026
10 Pack.			19.026-P10	19.026-P10
TIN		19.670	19.626	19.626
Recommended torque		35 Ncm	35 Ncm	35 Ncm
On abutment			19.009	19.009
10 Pack.			19.009-P10	19.009-P10
TIN			19.609	19.609
Recommended torque			15 Ncm	15 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.070	61.009	61.045
Ø 14mm	Ti		62.009	62.045
Ø 14mm	CoCr	71.070	71.009	71.045




CoCr BASE

		NNC	RN	WN
		3.5	4.1/4.8	4.8/6.5
			30.009	30.045
			31.009	31.045
	For welding techniques		30.109	30.145
			31.109	31.145
	On abutment		30.010	30.048
			31.010	31.048
For welding techniques			30.110	



CASTABLE

			28.009	28.045
	10 Pack.		11.009-P10	
			29.009	29.045
			12.009-P10	
	On abutment		28.010	28.046
	10 Pack.		11.010-P10	
			29.010	29.046
	On abutment			
	10 Pack.		12.010-P10	



CONICAL BL


compatible with
Straumann® Bone level

○ non-engaging
◻ engaging



HEALING ABUTMENT

	NP	RP	NP/RP
	NC/3.3	RC/4.1-4.8	NC/RC
GH 3,0mm	20.041/3	20.042/3	
GH 4,0mm	20.041/4	20.042/4	
GH 5,0mm	20.041/5	20.042/5	



SCAN ABUTMENTS



On abutment SRA

50.043	50.044	
		50.098



On abutment SRA

52.043	52.044	
		52.098



10 Pack.	○	53.100-P10	53.100-P10	
10 Pack.	◻	53.200-P10	53.200-P10	



TRANSFER

17.046	17.047
--------	--------

Transfers include two screws (for open and closed tray).



Platform Colour Coded

DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

OF SRA

*Other retention screws available, see page 59

NP	RP	NP/RP
NC/3.3	RC/4.1-4.8	NC/RC

14.043/D	14.044/D	
		14.098/D



Platform Colour Coded

ANALOGUE

10 Pack.

OF SRA

14.043	14.044	
14.043-P10	14.044-P10	

14.098



TEMPORARY ABUTMENT

On abutment SRA

○	24.043	24.044	
◻	25.043	25.044	
○			24.098
◻			25.098



DESSLoc®

GH 1,0mm		88.042/1
GH 2,0mm	88.041/2	88.042/2
GH 3,0mm	88.041/3	88.042/3
GH 4,0mm	88.041/4	88.042/4



MULTI-UNIT
ABUTMENTS

				NP	RP	NP/RP
				NC/3.3	RC/4.1-4.8	NC/RC
GH 1,5mm	0°	⦿		40.043/1	40.044/1	
GH 2,5mm	0°	⦿		40.043/2	40.044/2	
GH 3,5mm	0°	⦿		40.043/3	40.044/3	
GH 4,5mm	0°	⦿			40.044/4	
GH 2,5mm	17°	⦿		41.743/2	41.744/2	
GH 3,5mm	17°	⦿		41.743/3	41.744/3	
GH 2,5mm	30°	⦿		43.043/3	43.044/3	
GH 3,5mm	30°	⦿		43.043/4	43.044/4	
				4.8 MULTI-UNIT PLATFORM		




TI-BASE

GH 1,2mm	⦿	15.043	15.044
GH 3,0mm	⦿	15.343	15.344
GH 1,2mm	⦿	16.043	16.044
GH 3,0mm	⦿	16.343	16.344
ON SRA	⦿		15.098
	⦿		16.098



TI-BASE CASTABLES

5 Pack.	⦿	32.001-P5	32.001-P5
5 Pack.	⦿	33.001-P5	33.001-P5
ON SRA	⦿		32.098-P5
5 Pack.	⦿		33.098-P5




AURUMBase®

		NP	RP	NP/RP
		NC/3.3	RC/4.1-4.8	NC/RC
	⦿	36.043	36.044	
	⦿	37.043	37.044	

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.100-P5	33.100-P5
5 Pack.	10°	33.101-P5	33.101-P5
5 Pack.	20°	33.102-P5	33.102-P5




C-Base®

GH 0,7mm	⦿	58.043	58.044
GH 3,0mm	⦿	58.043/3	58.044/3
GH 0,7mm	⦿	59.043	59.044
GH 3,0mm	⦿	59.043/3	59.044/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



ELLIPTIBase®

A corner	⦿	56.043
B flat	⦿	57.043



Torx® Screw
AURUMBase® & ELLIPTIBase®

(replacement screw)	19.443	19.443
Recommended torque	35 Ncm	35 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



		NP	RP	NP/RP
ABUTMENTS		NC/3.3	RC/4.1-4.8	NC/RC
	0°	13.043	13.044	
	15°	22.043	22.044	



SCREW TORX®				
		19.032	19.033	19.037
10 Pack.		19.032-P10	19.033-P10	19.037-P10
TIN		19.632	19.633	19.637
Recommended torque		35 Ncm	35 Ncm	35 Ncm
On abutment SRA				19.098
10 Pack.				19.098-P10
TIN				19.698
Recommended torque				15 Ncm



PRE-MILLED BLANK				
Ø 10mm	Ti	61.043	61.044	
Ø 14mm	Ti	62.043	62.044	
Ø 14mm	CoCr	71.043	71.044	



CoCr BASE		NP	RP	NP/RP
		NC/3.3	RC/4.1-4.8	NC/RC
	○	30.043	30.044	
	⦿	31.043	31.044	
On abutment SRA	○			30.098
	⦿			31.098




CASTABLE				
	○	28.043	28.044	
10 Pack.		11.043-P10	11.044-P10	
	⦿	29.043	29.044	
10 Pack.		12.043-P10	12.044-P10	



CONICAL BLX


compatible with
Straumann® BLX

- non-engaging
- engaging



HEALING ABUTMENT	RB	WB	RB/WB
	3.5	4.5	
GH 0,8mm ^{A/B}		20.463/1	
GH 1,5mm ^A	20.461/2	20.461/2	20.461/2
GH 1,5mm ^B	20.462/2	20.462/2	20.462/2
GH 1,5mm ^{A/B}		20.463/2	
GH 2,5mm ^A	20.461/3	20.461/3	20.461/3
GH 2,5mm ^B	20.462/3	20.462/3	20.462/3
GH 3,5mm ^A	20.461/4	20.461/4	20.461/4
GH 3,5mm ^B	20.462/4	20.462/4	20.462/4


^Afor single crown restoration ^Bfor bridges and structures restorations



SCAN ABUTMENTS			
			52.062
	10 Pack.		53.100-P10
	10 Pack.		53.200-P10



TRANSFER			
			17.062
Transfers include two screws (for open and closed tray).			




Platform Colour Coded Retention screw Hex 1,27 mm*
*Other retention screws available, see page 59

DIGITAL ANALOGUE	RB	WB	RB/WB
	3.5	4.5	
	14.062/D	14.063/D	



Platform Colour Coded

ANALOGUE			
	14.062	14.063	



MULTI-UNIT ABUTMENTS				
GH 1,5mm	0°			40.062/1
GH 2,5mm	0°			40.062/2
GH 3,5mm	0°			40.062/3
GH 2,5mm	17°			41.762/2
GH 3,5mm	17°			41.762/3
GH 3,5mm	30°			43.062/3
GH 4,5mm	30°			43.062/4


4.8
MULTI-UNIT PLATFORM



			RB	WB	RB/WB
			3.5	4.5	
TI-BASE					
	GH 0,8mm	○		15.063	
	GH 1,5mm	○	15.062/2	15.063/2	15.062/2
	GH 2,5mm	○	15.062/3		15.062/3
	GH 3,5mm	○	15.062/4		15.062/4
	GH 0,8mm	⬢		16.063	
	GH 1,5mm	⬢	16.062/2	16.063/2	16.062/2
	GH 2,5mm	⬢	16.062/3		16.062/3
	GH 3,5mm	⬢	16.062/4		16.062/4



TI-BASE CASTABLES					
5 Pack.	○	32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.	⬢	33.001-P5	33.001-P5	33.001-P5	33.001-P5



AURUMBase®					
	GH 1,5mm	○	36.062	36.063	36.062
	GH 1,5mm	⬢	37.062	37.063	37.062



PLASTIC SLEEVES FOR AURUMBase®					
5 Pack.	0°		33.100-P5	33.100-P5	33.100-P5
5 Pack.	10°		33.101-P5	33.101-P5	33.101-P5
5 Pack.	20°		33.102-P5	33.102-P5	33.102-P5



Torx® Screw AURUMBase®

	RB	WB	RB/WB
	3.5	4.5	
(replacement screw)	19.441	19.441	19.441
Recommended torque	35 Ncm	35 Ncm	35 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW TORX®			
	19.054	19.054	19.054
TIN	19.700	19.700	19.700
Recommended torque	35 Ncm	35 Ncm	35 Ncm



PRE-MILLED BLANK			
Ø 10mm	Ti		61.062
Ø 14mm	Ti	62.063	
Ø 14mm	CoCr	71.063	



○ non-engaging
⬡ engaging



SCAN
ABUTMENTS



10 Pack. ○
10 Pack. ⬡

GM

52.035
53.100-P10
53.200-P10



TRANSFER

17.035

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

Platform Colour Coded Retention screw Hex 1,27 mm*
*Other retention screws available, see page 59

14.035/D



ANALOGUE

Platform Colour Coded

14.035



MULTI-UNIT
ABUTMENTS

GH 1,5mm 0° ○
GH 2,5mm 0° ○
GH 3,5mm 0° ○
GH 2,5mm 17° ⬡
GH 3,5mm 17° ⬡
GH 3,5mm 30° ⬡
GH 4,5mm 30° ⬡

GM

40.035/1
40.035/2
40.035/3
41.735/2
41.735/3
43.035/3
43.035/4



TI-BASE

GH 0,8mm ○
GH 1,5mm ○
GH 2,5mm ○
GH 0,8mm ⬡
GH 1,5mm ⬡
GH 2,5mm ⬡


15.035
15.235
15.335
16.035
16.235
16.335



TI-BASE CASTABLES

5 Pack. ○
5 Pack. ⬡

32.001-P5
33.001-P5



AURUMBase®

GH 2,5mm	○	36.035/3
GH 2,5mm	⬡	37.035/3

PLASTIC SLEEVES FOR AURUMBase®

	5 Pack.	0°	33.100-P5
	5 Pack.	10°	33.101-P5
	5 Pack.	20°	33.102-P5

Torx® Screw AURUMBase®



(replacement screw)	19.444
Recommended torque	20 Ncm

Torx® BALL SCREWDRIVER



20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW NEO

	19.051
<i>Recommended torque</i>	<i>20 Ncm</i>




PRE-MILLED BLANK

Ø 10mm	Ti	61.035
Ø 14mm	Ti	62.035
Ø 14mm	CoCr	71.035

DESSLoc®

DESSLoc® is the best alternative for your removable overdenture projects. Its Zirconium Nitride surface coating improves the overall mechanical resistance whilst reducing gingival inflammation and plaque adhesion.




 non-engaging

 engaging



HEALING
ABUTMENT

	NP	RP	WP
	3.5	4.5	5.7
GH 3,0mm	20.017/3	20.018/3	20.019/3
GH 4,0mm	20.017/4	20.018/4	20.019/4
GH 5,0mm	20.017/5	20.018/5	20.019/5



SCAN
ABUTMENTS

		50.017	50.018	50.019
		52.017	52.018	52.019
	10 Pack. 	53.100-P10	53.100-P10	53.100-P10
	10 Pack. 	53.200-P10	53.200-P10	53.200-P10



TRANSFER

	17.017	17.018	17.019
Transfers include two screws (for open and closed tray).			



DIGITAL
ANALOGUE

	Platform Colour Coded	Retention screw Hex 1,27 mm*	14.017/D	14.018/D	14.019/D
*Other retention screws available, see page 59					



ANALOGUE

	Platform Colour Coded	10 Pack.	14.017	14.018	14.019
			14.017-P10	14.018-P10	



TEMPORARY
ABUTMENT

	NP	RP	WP
	3.5	4.5	5.7
	24.017	24.018	24.019
	25.017	25.018	25.019



DESSLoc®

GH 1,0mm	88.017/1	88.018/1	88.019/1
GH 2,5mm	88.017/2	88.018/2	88.019/2
GH 3,5mm	88.017/3	88.018/3	88.019/3
GH 4,5mm	88.017/4	88.018/4	88.019/4
GH 5,5mm	88.017/5	88.018/5	




MULTI-UNIT
ABUTMENTS

GH 1,0mm	0°		40.019/1	40.020/1
GH 2,0mm	0°		40.019/2	40.020/2
GH 3,0mm	0°		40.019/3	40.020/3
GH 4,0mm	0°		40.019/4	
GH 5,0mm	0°		40.019/5	
GH 2,5mm	17°		41.719/2	41.720/2
GH 3,5mm	17°		41.719/3	41.720/3
GH 2,5mm	30°		43.019/3	43.020/3
GH 3,5mm	30°		43.019/4	43.020/4


4.8
MULTI-UNIT PLATFORM



			NP	RP	WP
			3.5	4.5	5.7
TI-BASE					
	GH 0,3mm	○		15.018	15.019
	GH 0,8mm	○	15.017		
	GH 1,5mm	○	15.217	15.218	15.219
	GH 3,0mm	○	15.317	15.318	15.319
	GH 0,3mm	○		16.018	16.019
	GH 0,8mm	○	16.017		
	GH 1,5mm	○	16.217	16.218	16.219
	GH 3,0mm	○	16.317	16.318	16.319
TI-BASE CASTABLES					
	5 Pack.	○	32.001-P5	32.001-P5	32.001-P5
	5 Pack.	○	33.001-P5	33.001-P5	33.001-P5




AURUMBase®					
		○	36.017	36.018	36.019
		○	37.017	37.018	37.019
PLASTIC SLEEVES FOR AURUMBase®					
	5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5



			NP	RP	WP
			3.5	4.5	5.7
C-Base®					
	GH 0,5mm	○	58.017	58.018	58.019
	GH 1,5mm	○	58.017/2	58.018/2	58.019/2
	GH 3,0mm	○	58.017/3	58.018/3	58.019/3
	GH 0,5mm	○	59.017	59.018	59.019
	GH 1,5mm	○	59.017/2	59.018/2	59.019/2
	GH 3,0mm	○	59.017/3	59.018/3	59.019/2

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



ELLIPTIBase®			
A corner	◇		56.017
B flat	□		57.017



Torx® Screw AURUMBase® & ELLIPTIBase®					
(replacement screw)			19.417	19.417	19.417
Recommended torque			30 Ncm	30 Ncm	30 Ncm



Torx® BALL SCREWDRIVER		20mm	24mm	30mm	35mm
		DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS					
	0°		13.017	13.018	13.019
	15°		22.017	22.018	22.019
	25°		23.017	23.018	



SCREW HEX 1,27mm

	NP	RP	WP
	3.5	4.5	5.7
	19.012	19.012	19.012
10 Pack.	19.012-P10	19.012-P10	19.012-P10
Recommended torque	30 Ncm	30 Ncm	30 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.017	61.018	61.019
Ø 14mm	Ti	62.017	62.018	62.019
Ø 14mm	CoCr	71.017	71.018	71.019



CoCr BASE

	◯	30.017	30.018	30.019
	◯	31.017	31.018	31.019



CASTABLE

	◯	28.017	28.018	28.019
10 Pack.		11.017-P10	11.018-P10	
	◯	29.017	29.018	29.019
10 Pack.		12.017-P10	12.018-P10	


Digital analogue

Our analogues have been carefully designed to match the implant connection to perfection. Manufactured either in stainless steel AISI 303 or anodized titanium grade 5 ELI.



- Color coded to match implant platform
- Precise connection guaranteed
- For 3D and traditional clay models
- Digital analogue includes HEX 1,20 screw

Retention screws	
HEX 1,00	19.707-P5
HEX 1,18	19.708-P5
HEX 1,20	19.709-P5
HEX 1,22	19.710-P5
HEX 1,25	19.711-P5
HEX 1,27	19.712-P5
HEX 1,28	19.713-P5
HEX 1,40	19.714-P5
Unigrip®	19.715-P5
Torx®	19.716-P5
NE	19.717-P5



CONIC EVO

compatible with
Astra Tech implant system™ EV

- non-engaging
- ⊗ engaging



HEALING ABUTMENT

	EV/3.0	EV/3.6	EV/4.2	EV/4.8	EV/5.4
	3.0	3.6	4.2	4.8	5.4
GH 3,5mm	20.406/3	20.059/3	20.060/3	20.061/3	20.407/3
GH 4,5mm	20.406/4	20.059/4	20.060/4	20.061/4	20.407/4
GH 6,5mm	20.406/6	20.059/6	20.060/6	20.061/6	20.407/6



SCAN ABUTMENTS

		50.059	50.060	50.061
		52.059	52.060	52.061
	10 / 5 Pack.	53.100-P10	53.110-P5	53.148-P5
	10 / 5 Pack.	53.200-P10	53.210-P5	53.248-P5
On uniabutment®	○	50.058		
On uniabutment®	⊗	52.058		



TRANSFER

		17.059	17.060	17.061
On uniabutment®	17.058			

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE

	EV/3.0	EV/3.6	EV/4.2	EV/4.8	EV/5.4
	3.0	3.6	4.2	4.8	5.4
Retention screw Hex 1,27 mm*		14.059/D	14.060/D	14.061/D	
On uniabutment®		14.058/D			

*Other retention screws available, see page 59



ANALOGUE

		14.059	14.060	14.061
On uniabutment®	14.058			



TEMPORARY ABUTMENT

	○	24.059	24.060	24.061
	⊗	25.059	25.060	25.061
On uniabutment®	○	24.058		



DESSLoc®

GH 1,0mm	88.059/1	88.060/1	88.061/1
GH 2,0mm	88.059/2	88.060/2	88.061/2
GH 3,0mm	88.059/3	88.060/3	88.061/3
GH 4,0mm	88.059/4	88.060/4	88.061/4



UNIABUTMENT®

	EV/3.0	EV/3.6	EV/4.2	EV/4.8	EV/5.4
	3.0	3.6	4.2	4.8	5.4
GH 1,0mm		40.059/1	40.060/1	40.061/1	
GH 2,0mm		40.059/2	40.060/2	40.061/2	
GH 3,0mm		40.059/3	40.060/3	40.061/3	
GH 5,0mm		40.059/5	40.060/5	40.061/5	

UniAbutment® needs a 1,75mm hex screwdriver



MULTI-UNIT
ABUTMENTS

GH 1,5mm	0°		40.027/1	40.036/1	40.037/1
GH 2,5mm	0°		40.027/2	40.036/2	40.037/2
GH 3,5mm	0°		40.027/3	40.036/3	40.037/3
GH 3,5mm	17°		41.727/3	41.736/3	41.737/3
GH 4,5mm	17°		41.727/4	41.736/4	41.737/4
GH 4,5mm	30°		43.027/4	43.036/4	43.037/4
GH 5,5mm	30°		43.027/5	43.036/5	43.037/5

4.8
MULTI-UNIT PLATFORM



TI-BASE


GH 1,0mm		15.059	15.060	15.061	15.407
GH 3,0mm		15.359	15.360	15.361	
GH 1,0mm		16.059	16.060	16.061	16.407
GH 3,0mm		16.359	16.360	16.361	

On uniabutment®  15.058



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5
5 Pack.		32.001-P5			



AURUMBase®

	EV/3.0	EV/3.6	EV/4.2	EV/4.8	EV/5.4
	3.0	3.6	4.2	4.8	5.4
		36.059	36.060	36.061	
		37.059	37.060	37.061	
On uniabutment® 	36.058				

PLASTIC SLEEVES FOR AURUMBase®

5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5




C-Base®

GH 1,0mm		58.059	58.060	58.061
GH 3,0mm		58.059/3	58.060/3	58.061/3
GH 1,0mm		59.059	59.060	59.061
GH 3,0mm		59.059/3	59.060/3	59.061/3


Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



ELLIPTIBase®

A corner		56.406
B flat		57.406



Torx® Screw
AURUMBase® & ELLIPTIBase®

(replacement screw)	19.441	19.460	19.461
Recommended torque	25 Ncm	25 Ncm	25 Ncm
On uniabutment®	19.331		
Recommended torque	15 Ncm		



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW HEX 1,27mm

		EV/3.0	EV/3.6	EV/4.2	EV/4.8	EV/5.4
TIN		3.0	3.6	4.2	4.8	5.4
			19.059	19.060	19.061	
			19.659	19.660	19.661	
Recommended torque			25 Ncm	25 Ncm	24 Ncm	
On uniabutment®	19.058					
TIN	19.658					
Recommended torque			25 Ncm	25 Ncm	24 Ncm	



PRE-MILLED BLANK

Ø 10mm	Ti	61.406	61.059	61.060	61.061	
Ø 14mm	Ti		62.059	62.060	62.061	62.407
Ø 14mm	CoCr		71.059	71.060	71.061	71.407



CoCr BASE

		30.059	30.060	30.061	
		31.059	31.060	31.061	
On uniabutment®		30.058			




CASTABLE

		28.059	28.060	28.061	
On uniabutment®		28.058			

DESS® Multi-Unit

With 14 different connections available, the standardization with a unique prosthetic platform can be achieved allowing you to have a full-arch rehab with up to 4 different implant compatibilities.



INTERNAL HEX CONIC

compatible with
Astra Tech Osseospeed™

○ non-engaging
○ engaging




HEALING ABUTMENT

			NP	RP	WP
			3.0	3.5-4.0	4.5-5.0
	GH 3,0mm			20.024/3	20.025/3
	GH 4,0mm			20.024/4	20.025/4
	GH 5,0mm			20.024/5	20.025/5
Anatomic	GH 4,0mm	Ø 4,5		20.424/4	
	GH 4,0mm	Ø 5,5		20.524/4	20.525/4
	GH 4,0mm	Ø 6,5			20.625/4



SCAN ABUTMENTS

			50.023	50.024	50.025
			52.023	52.024	52.025
	10 Pack.	○	53.100-P10	53.100-P10	53.100-P10
	10 Pack.	○	53.200-P10	53.200-P10	53.200-P10
On uniabutment®	20°		50.036		
	45°		50.037		
On uniabutment®	20°		52.036		
	45°		52.037		



TRANSFER

			17.023	17.024	17.025
On uniabutment®	20°		17.036		
	45°		17.037		



Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE



Platform Colour Coded

			NP	RP	WP
			3.0	3.5-4.0	4.5-5.0
Retention screw Hex 1,27 mm*			14.023/D	14.024/D	14.025/D
Of uniabutment®	20°		14.036/D		
	45°		14.037/D		

*Other retention screws available, see page 59




ANALOGUE



Platform Colour Coded

			14.023	14.024	14.025
			10 Pack.	14.024-P10	14.025-P10
Of uniabutment®	20°		14.036		
	10 Pack.		14.036-P10		
	45°		14.037		
	10 Pack.		14.037-P10		



TEMPORARY ABUTMENT

			○	24.024	24.025
			○	25.024	25.025
On uniabutment®	20°	○	24.028		
	45°	○	24.029		



DESSLoc®

		NP	RP	WP
		3.0	3.5-4.0	4.5-5.0
GH 1,0mm		88.023/1	88.024/1	88.025/1
GH 2,0mm		88.023/2	88.024/2	88.025/2
GH 3,0mm		88.023/3	88.024/3	88.025/3
GH 4,0mm		88.023/4	88.024/4	88.025/4


UNIABUTMENT®

GH 0,5mm	20°		40.013/0	40.015/0
GH 1,0mm	20°		40.013/1	40.015/1
GH 2,0mm	20°		40.013/2	40.015/2
GH 4,0mm	20°		40.013/4	40.015/4
GH 6,0mm	20°		40.013/6	40.015/6
GH 0,5mm	45°		40.014/0	40.016/0
GH 1,0mm	45°		40.014/1	40.016/1
GH 2,0mm	45°		40.014/2	40.016/2
GH 4,0mm	45°		40.014/4	40.016/4
GH 6,0mm	45°		40.014/6	40.016/6

MULTI-UNIT
ABUTMENTS

GH 1,5mm	0°	○	40.024/1	40.025/1
GH 2,5mm	0°	○	40.024/2	40.025/2
GH 3,5mm	0°	○	40.024/3	40.025/3
GH 3,5mm	17°	○	41.724/3	41.725/3
GH 4,5mm	17°	○	41.724/4	41.725/4
GH 4,5mm	30°	○	43.024/4	43.025/4
GH 5,5mm	30°	○	43.024/5	43.025/5

4.8
MULTI-UNIT PLATFORM



TI-BASE

		NP	RP	WP
		3.0	3.5-4.0	4.5-5.0
GH 0,9mm	○			15.025
GH 1,0mm	○	15.023		
GH 1,1mm	○		15.024	
GH 2,0mm	○	15.023/2	15.024/2	15.025/2
GH 3,0mm	○	15.023/3	15.324	15.325
GH 0,9mm	○			16.025
GH 1,0mm	○	16.023		
GH 1,1mm	○		16.024	
GH 2,0mm	○	16.023/2	16.024/2	16.025/2
GH 3,0mm	○	16.023/3	16.324	16.325
On uniabutment®	20°	○	15.036	
	45°	○	15.037	

TI-BASE CASTABLES



5 Pack.	○	32.001-P5	32.001-P5	32.001-P5
5 Pack.	○	33.001-P5	33.001-P5	33.001-P5
On uniabutment®		32.001-P5		


AURUMBase®

	○		36.024	36.025
	○		37.024	37.025

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°		33.100-P5	33.100-P5
5 Pack.	10°		33.101-P5	33.101-P5
5 Pack.	20°		33.102-P5	33.102-P5



		NP	RP	WP
		3.0	3.5-4.0	4.5-5.0
C-Base®				
	GH 0,8mm		58.024	
	GH 0,9mm			58.025
	GH 3,0mm		58.024/3	58.025/3
	GH 0,8mm		59.024	
	GH 0,9mm			59.025
	GH 3,0mm		59.024/3	59.025/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



Torx® Screw AURUMBase®			
(replacement screw)		19.441	19.425
Recommended torque		20 Ncm	25 Ncm



Torx® BALL SCREWDRIVER		20mm	24mm	30mm	35mm
		DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS			
	0°	13.024	13.025
	15°	22.024	22.025
	25°	23.024	



		NP	RP	WP
		3.0	3.5-4.0	4.5-5.0
SCREW HEX 1,27mm				
On uniabutment®		19.015	19.016	19.017
10 Pack.			19.016-P10	19.017-P10
TIN		19.615	19.616	19.617
Recommended torque		15 Ncm	20 Ncm	25 Ncm
		19.018		
10 Pack.			19.018-P10	
TIN			19.618	
Recommended torque			15 Ncm	

PRE-MILLED BLANK				
Ø 10mm	Ti	61.023	61.024	61.025
Ø 14mm	Ti		62.024	62.025
Ø 14mm	CoCr		71.024	71.025

CoCr BASE			
			
		30.024	30.025
		31.024	31.025
On uniabutment®		30.036	
		30.037	

CASTABLE				
			28.024	28.025
10 Pack.			11.024-P10	11.025-P10
			29.024	29.025
10 Pack.			12.024-P10	12.025-P10
On uniabutment®	20°		28.036	
10 Pack.			11.036-P10	
On uniabutment®	45°		28.037	
10 Pack.			11.037-P10	



INTERNAL HEX FD

compatible with
Xive®


 non-engaging

 engaging



HEALING ABUTMENT

	NP	RP	WP	5.5
	3.4	3.8	4.5	5.5
GH 3,0mm	20.036/3	20.037/3	20.038/3	
GH 4,0mm	20.036/4	20.037/4	20.038/4	
GH 5,0mm	20.036/5	20.037/5	20.038/5	




SCAN ABUTMENTS

		50.038	50.039	50.040	50.071
		52.038	52.039	52.040	52.071
	10 Pack.	53.100-P10	53.100-P10	53.100-P10	53.100-P10
	10 Pack.	53.200-P10	53.200-P10	53.200-P10	53.200-P10



TRANSFER

17.038	17.039	17.040	17.071
Transfers include two screws (for open and closed tray).			



Platform Colour Coded

DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

NP	RP	WP	5.5
3.4	3.8	4.5	5.5
14.038/D	14.039/D	14.040/D	14.071/D



Platform Colour Coded

ANALOGUE

14.038	14.039	14.040	14.071
10 Pack.	14.039-P10		



TEMPORARY ABUTMENT

	24.038	24.039	24.040
	25.038	25.039	25.040



DESSLoc®

GH 1,0mm	88.036/1	88.037/1
GH 2,0mm	88.036/2	88.037/2
GH 3,0mm	88.036/3	88.037/3
GH 4,0mm	88.036/4	88.037/4




TI-BASE

			NP	RP	WP	5.5
			3.4	3.8	4.5	5.5
GH 0,3mm	○		15.038	15.039	15.040	15.071
GH 1,5mm	○		15.238	15.239	15.240	15.271
GH 3,0mm	○		15.338	15.339	15.340	15.371
GH 0,3mm	⬡		16.038	16.039	16.040	16.071
GH 1,5mm	⬡		16.238	16.239	16.240	16.271
GH 3,0mm	⬡		16.338	16.339	16.340	16.371



TI-BASE CASTABLES

5 Pack.	○	32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.	⬡	33.001-P5	33.001-P5	33.001-P5	33.001-P5




AURUMBase®

○	36.038	36.039	36.040
⬡	37.038	37.039	37.040

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
5 Pack.	10°	33.100-P5	33.101-P5	33.101-P5
5 Pack.	20°	33.100-P5	33.102-P5	33.102-P5



C-Base®

			NP	RP	WP	5.5
			3.4	3.8	4.5	5.5
GH 0,3mm	○		58.038	58.039		
GH 0,6mm	○				58.040	58.071
GH 1,5mm	○		58.038/2	58.039/2	58.040/2	58.071/2
GH 3,0mm	○		58.038/3	58.039/3	58.040/3	58.071/3
GH 0,3mm	⬡		59.038	59.039		
GH 0,6mm	⬡				59.040	59.071
GH 1,5mm	⬡		59.038/2	59.039/2	59.040/2	59.071/2
GH 3,0mm	⬡		59.038/3	59.039/3	59.040/3	59.071/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



Torx® Screw AURUMBase®

(replacement screw)	19.438	19.438	19.438
Recommended torque	24 Ncm	24 Ncm	24 Ncm



Torx® BALL SCREWDRIVER

	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS

0°	13.038	13.039	13.040
15°	22.038	22.039	22.040
25°		23.039	



SCREW HEX 1,22mm

	NP	RP	WP	5.5
	3.4	3.8	4.5	5.5
	19.024	19.024	19.024	19.024
10 Pack.	19.024-P10	19.024-P10	19.024-P10	19.024-P10
TIN	19.624	19.624	19.624	19.624
Recommended torque	24 Ncm	24 Ncm	24 Ncm	24 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.038	61.039	61.040	61.071
Ø 14mm	Ti		62.039	62.040	62.071
Ø 14mm	CoCr	71.038	71.039	71.040	71.071



CoCr BASE

○	30.038	30.039	30.040
⬡	31.038	31.039	31.040




CASTABLE

	○	28.038	28.039	28.040
10 Pack.			11.039-P10	
		29.038	29.039	29.040
10 Pack.	⬡		12.039-P10	


We are DESS

DESS® offers a variety of dental prosthetic solutions for nearly all the connections available in the market. DESS® product quality is as good as or better than the quality of well-known brands. The “pure switch concept” on which our products are conceived assures full compatibility with the implant brand chosen.



Reasons to switch to DESS®

- **One Stop Shop**
31 compatibilities with the most popular implant connections
- **European Quality**
ISO 9001, ISO 13485, CE, FDA, Health Canada
- **Lifetime Warranty**
All DESS® manufactured products have a lifetime warranty*



INTERNAL ANK

compatible with
Dentsply Ankylos® C/X

- non-engaging
- ⬡ engaging



HEALING ABUTMENT

GH 1,5mm	20.046/1
GH 2,5mm	20.046/3



SCAN ABUTMENTS

	○	50.049
	⬡	52.049
	10 Pack.	53.100-P10
	10 Pack.	53.200-P10



TRANSFER

GH 1,5mm	17.049
GH 2,5mm	17.250

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE

Platform Colour Coded	Retention screw Hex 1,27 mm*	Ø 3,5	14.048/D
		Ø 4,5	14.049/D

*Other retention screws available, see page 59



ANALOGUE

Platform Colour Coded	Ø 3,5	14.048
	10 Pack.	14.048-P10
	Ø 4,5	14.049
	10 Pack.	14.049-P10



TEMPORARY ABUTMENT

GH 1,5mm	○	24.048
GH 2,5mm	○	24.250
GH 1,5mm	⬡	25.048
GH 2,5mm	⬡	25.250



DESSLoc®

GH 2,0mm	88.046/2
GH 3,0mm	88.046/3
GH 4,0mm	88.046/4



TI-BASE

GH 1,5mm		15.050
GH 2,5mm		15.250
GH 1,5mm		16.050
GH 2,5mm		16.250



TI-BASE CASTABLES

5 Pack.		32.001-P5
5 Pack.		33.001-P5



C-Base®

GH 1,4mm		58.050
GH 2,5mm		58.050/2
GH 1,4mm		59.050
GH 2,5mm		59.050/2

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



Replacement Screw for C-Base®

	19.702
<i>Recommended torque</i>	15 Ncm



ABUTMENTS

GH 1,5mm	0°	13.047
GH 2,5mm	0°	13.250
GH 1,5mm	15°	22.048
GH 2,5mm	15°	22.250
GH 1,5mm	25°	23.048
GH 2,5mm	25°	23.250



SCREW HEX 1,00mm

GH 1,5mm	19.048
10 Pack.	19.048-P10
Tin	19.648
GH 2,5mm	19.050

Recommended torque

15 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.250
Ø 14mm	CoCr	71.250




CoCr BASE

GH 1,5mm		30.050
GH 2,5mm		30.250
GH 1,5mm		31.050
GH 2,5mm		31.250



CASTABLE


GH 1,5mm		28.050
GH 2,5mm		28.250



CONIC OSS

compatible with
Osstem® TS / HIOSSEN® ET

- non-engaging
- engaging




SCAN
ABUTMENTS



- 10 Pack.
- 10 Pack.

	Mini	STD
	Mini	Regular
	50.072	50.073
	52.072	52.073
	53.100-P10	53.100-P10
	53.200-P10	53.200-P10



TRANSFER

17.072	17.073
--------	--------

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

Platform Colour Coded

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59


14.072/D	14.073/D
----------	----------



ANALOGUE

Platform Colour Coded

14.072	14.073
--------	--------



TEMPORARY
ABUTMENT

-
-

	Mini	STD
	Mini	Regular
	24.072	24.073
	25.072	25.073



MULTI-UNIT
ABUTMENTS

GH 1,5mm	0°	
GH 2,5mm	0°	
GH 3,5mm	0°	
GH 2,5mm	17°	
GH 3,5mm	17°	
GH 3,5mm	30°	
GH 4,5mm	30°	

4.8
MULTI-UNIT PLATFORM



TI-BASE


GH 1,0mm	
GH 3,0mm	
GH 1,0mm	
GH 3,0mm	



TI-BASE CASTABLES


- 5 Pack.
- 5 Pack.

32.001-P5	32.001-P5
33.001-P5	33.001-P5



AURUMBase®

			36.072	36.073
			37.072	37.073
PLASTIC SLEEVES FOR AURUMBase®				
	5 Pack.	0°	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5



C-Base®

GH 1,0mm		58.072	58.073
GH 3,0mm		58.072/3	58.073/3
GH 1,0mm		59.072	59.073
GH 3,0mm		59.072/3	59.073/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.



Torx® Screw AURUMBase®

(replacement screw)	19.444	19.461
Recommended torque	20 Ncm	25 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW HEX 1,22mm

		Mini	STD
		Mini	Regular
		19.072	19.073
	DLC	19.272	19.273
	TIN	19.672	19.673
Recommended torque		20 Ncm	30 Ncm




PRE-MILLED BLANK

Ø 10mm	Ti	61.072	61.073
Ø 14mm	Ti	62.072	62.073
Ø 14mm	CoCr	71.072	71.073



CoCr BASE

		30.072	30.073
		31.072	31.073



INTERNAL CAM


compatible with
Camlog®

-  non-engaging
-  engaging



HEALING ABUTMENT

	3.3	3.8	4.3	5.0
	3.3	3.8	4.3	5.0
GH 3,0mm		20.065/3	20.066/3	
GH 4,0mm		20.065/4	20.066/4	
GH 5,0mm		20.065/5	20.066/5	



SCAN ABUTMENTS



	50.064	50.065	50.066	50.067
	52.064	52.065	52.066	52.067
10 Pack. 	53.100-P10	53.100-P10	53.100-P10	53.100-P10
10 Pack. 	53.200-P10	53.200-P10	53.200-P10	53.200-P10



TRANSFER

17.064	17.065	17.066	17.067
Transfers include two screws (for open and closed tray).			



Platform Colour Coded

DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

3.3	3.8	4.3	5.0
3.3	3.8	4.3	5.0
14.064/D	14.065/D	14.066/D	14.067/D



Platform Colour Coded

ANALOGUE

14.064	14.065	14.066	14.067
--------	--------	--------	--------



TEMPORARY ABUTMENT

	24.065	24.066
	25.065	25.066



DESSLoc®


GH 1,0mm	88.065/1	88.066/1
GH 2,0mm	88.065/2	88.066/2
GH 3,0mm	88.065/3	88.066/3
GH 4,0mm	88.065/4	88.066/4



TI-BASE			3.3	3.8	4.3	5.0
			3.3	3.8	4.3	5.0
	GH 0,3mm			15.065	15.066	15.067
	GH 0,5mm		15.064			
	GH 1,5mm			15.265	15.266	
	GH 3,0mm			15.365	15.366	
	GH 0,3mm			16.065	16.066	16.067
	GH 0,5mm		16.064			
	GH 1,5mm			16.265	16.266	
	GH 3,0mm			16.365	16.366	



TI-BASE CASTABLES						
5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5	
5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5	



AURUMBase®					
			36.065	36.066	
			37.065	37.066	



PLASTIC SLEEVES FOR AURUMBase®					
5 Pack.	0°		33.100-P5	33.100-P5	
5 Pack.	10°		33.101-P5	33.101-P5	
5 Pack.	20°		33.102-P5	33.102-P5	



Torx® Screw AURUMBase®				
(replacement screw)			19.444	19.444
Recommended torque			20 Ncm	20 Ncm



Torx® BALL SCREWDRIVER	3.3	3.8	4.3	5.0
	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW HEX 1,27mm					
		19.044	19.044	19.044	19.045
	TIN	19.644	19.644	19.644	19.645
Recommended torque		20 Ncm	20 Ncm	20 Ncm	20 Ncm




PRE-MILLED BLANK					
Ø 10mm	Ti	61.064	61.065	61.066	61.067
Ø 14mm	Ti		62.065	62.066	62.067
Ø 14mm	CoCr	71.064	71.065	71.066	71.067



CoCr BASE					
			30.065	30.066	
			31.065	31.066	



CASTABLE					
			28.065	28.066	



INTERNAL CON

compatible with
Conelog®

-  non-engaging
-  engaging



SCAN ABUTMENTS



- 10 Pack. 
- 10 Pack. 

	NP	RP	WP
	3.3	3.8/4.3	5.0
	52.074	52.075	52.076
	53.100-P10	53.100-P10	53.100-P10
	53.200-P10	53.200-P10	53.200-P10



TRANSFER

17.074	17.075	17.076
--------	--------	--------

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

14.074/D	14.075/D	14.076/D
----------	----------	----------



ANALOGUE



Platform Colour Coded

14.074	14.075	14.076
--------	--------	--------




TI-BASE

		NP	RP	WP
		3.3	3.8/4.3	5.0
GH 1,4mm		15.074	15.075	15.076
GH 3,0mm		15.074/3	15.075/3	15.076/3
GH 1,4mm		16.074	16.075	16.076
GH 3,0mm		16.074/3	16.075/3	16.076/3



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5




AURUMBase®

	36.074	36.075	36.076
	37.074	37.075	37.076

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5



C-Base®

GH 1,0mm	
GH 2,1mm	
GH 1,0mm	
GH 2,1mm	

Includes screw, same as standard screw.

NP	RP	WP
3.3	3.8/4.3	5.0
58.074	58.075	58.076
58.074/2	58.075/2	58.076/2
59.074	59.075	59.076
59.074/2	59.075/2	59.076/2

Libraries available for EXOCAD & 3Shape.



Replacement Screw for C-Base®

19.052	19.052	19.053
20 Ncm	20 Ncm	20 Ncm

Torx® Screw AURUMBase®



(replacement screw)	19.474	19.474	19.476
Recommended torque	20 Ncm	20 Ncm	20 Ncm

Torx® BALL SCREWDRIVER



20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW HEX 1,27mm

TIN


Recommended torque

NP	RP	WP
3.3	3.8/4.3	5.0
19.074	19.074	19.076
19.674	19.674	19.676
20 Ncm	20 Ncm	20 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.074	61.075	61.076
Ø 14mm	Ti	62.074	62.075	62.076
Ø 14mm	CoCr	71.074	71.075	71.076



INTERNAL HEX BH


compatible with
BioHorizons® internal

 non-engaging

 engaging



HEALING ABUTMENT	3.0	NP	RP	WP
	3.0	3.5	4.5	5.7
GH 3,0mm		20.087/3	20.088/3	20.089/3
GH 4,0mm		20.087/4	20.088/4	20.089/4
GH 5,0mm		20.087/5	20.088/5	20.089/5



SCAN ABUTMENTS					
		52.099	52.087	52.088	52.089
	10 Pack. 		53.100-P10	53.100-P10	53.100-P10
	10 Pack. 		53.200-P10	53.200-P10	53.200-P10



TRANSFER				
		17.087	17.088	17.089
		Transfers include two screws (for open and closed tray).		



	Platform Colour Coded	Retention screw Hex 1,27 mm*	14.099/D	14.087/D	14.088/D	14.089/D
		*Other retention screws available, see page 59				

*Other retention screws available, see page 59



ANALOGUE	3.0	NP	RP	WP
	3.0	3.5	4.5	5.7
	14.099	14.087	14.088	14.089





DESSLoc®			
GH 1,0mm		88.087/1	88.088/1
GH 2,5mm		88.087/2	88.088/2
GH 3,5mm		88.087/3	88.088/3
GH 4,5mm		88.087/4	88.088/4
GH 5,5mm		88.087/5	



TI-BASE				
GH 0,3mm			15.088	15.089
GH 0,8mm		15.087		
GH 0,3mm			16.088	16.089
GH 0,8mm		16.087		



TI-BASE CASTABLES				
5 Pack.		32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5



MULTI-UNIT
ABUTMENTS

			3.0	NP	RP	WP
			3.0	3.5	4.5	5.7
GH 1,0mm	0°	◯		40.087/1	40.088/1	
GH 2,0mm	0°	◯		40.087/2	40.088/2	
GH 3,0mm	0°	◯		40.087/3	40.088/3	
GH 4,0mm	0°	◯		40.087/4		
GH 5,0mm	0°	◯		40.087/5		
GH 2,5mm	17°	◻		41.787/2	41.720/2	
GH 3,5mm	17°	◻		41.787/3	41.720/3	
GH 3,5mm	30°	◻		43.087/3	43.020/3	
GH 4,5mm	30°	◻		43.087/4	43.020/4	

4.8
MULTI-UNIT PLATFORM

AURUMBase®


	◯	36.087	36.088	36.089
	◻	37.087	37.088	37.089

PLASTIC SLEEVES FOR AURUMBase®

	5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5

ELLIPTIBase®

A corner	◻	56.099
B flat	◻	57.099



C-Base®

		3.0	NP	RP	WP
		3.0	3.5	4.5	5.7
GH 0,5mm	◯		58.087	58.088	58.089
GH 3,0mm	◯		58.087/3	58.088/3	58.089/3
GH 0,5mm	◻		59.087	59.088	59.089
GH 3,0mm	◻		59.087/3	59.088/3	59.089/3

Includes screw, same as standard screw.

Libraries available for EXOCAD & 3Shape.

Torx® Screw
AURUMBase® & ELLIPTIBase®

(replacement screw)	19.400	19.417	19.417	19.417
Recommended torque	30 Ncm	30 Ncm	30 Ncm	30 Ncm

Torx® BALL SCREWDRIVER

	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06

SCREW HEX 1,27mm

	19.043	19.043	19.043	19.043
TIN	19.643	19.643	19.643	19.643
Recommended torque	30 Ncm	30 Ncm	30 Ncm	30 Ncm

PRE-MILLED BLANK

Ø 10mm	Ti	61.099	61.087	61.088	61.089
Ø 14mm	Ti			62.088	62.089

EXTERNAL
HEX BH

compatible with
BioHorizons® External


 non-engaging

 engaging



HEALING
ABUTMENT

	NP	RP	WP
	3.7	4.2	5.2
GH 3,0mm	20.020/3	20.021/3	20.022/3
GH 4,0mm	20.020/4	20.021/4	20.022/4
GH 5,0mm	20.020/5	20.021/5	20.022/5



SCAN
ABUTMENTS




		50.020	50.021	50.022
		52.020	52.021	52.022
10 Pack.		53.100-P10	53.100-P10	53.100-P10
10 Pack.		53.200-P10	53.200-P10	53.200-P10



TRANSFER

17.020	17.021	17.022
--------	--------	--------

Transfers include two screws (for open and closed tray).



Platform Colour Coded

DIGITAL
ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

NP	RP	WP
3.7	4.2	5.2
14.020/D	14.021/D	14.022/D



Platform Colour Coded

ANALOGUE

10 Pack.

14.020	14.021	14.022
	14.021-P10	



TEMPORARY
ABUTMENT



24.020	24.021	24.022
25.020	25.021	25.022



TI-BASE

GH 0,3mm		15.051	15.052	15.053
GH 0,3mm		16.051	16.052	16.053



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5



ABUTMENTS

		NP	RP	WP
		3.7	4.2	5.2
	0°	13.020	13.021	13.022
	15°	22.020	22.021	



SCREW HEX 1,27mm

		19.013	19.014	19.014
10 Pack.			19.014-P10	19.014-P10
TIN		19.613	19.614	19.614
Recommended torque		30 Ncm	30 Ncm	30 Ncm




PRE-MILLED BLANK


Ø 10mm	Ti	61.020	61.021	61.022
Ø 14mm	Ti	62.020	62.021	62.022
Ø 14mm	CoCr	71.020	71.021	71.022



CASTABLE

10 Pack.		28.020	28.021	28.022
			11.021-P10	
10 Pack.		29.020	29.021	29.022
			12.021-P10	

 C-Base





INTERNAL HEX MI


compatible with
Mis® Seven

 non-engaging

 engaging



HEALING ABUTMENT	NP	SP	WP	Multi-Unit
	3.3	3.75-4.2	5.0-6.0	4.8
GH 3,0mm		20.017/3	20.018/3	
GH 4,0mm		20.017/4	20.018/4	
GH 5,0mm		20.017/5	20.018/5	





SCAN ABUTMENTS

		50.056	50.017	50.018	
		52.056	52.017	52.018	52.085
	10 Pack. 	53.100-P10	53.100-P10	53.100-P10	
	10 Pack. 	53.200-P10	53.200-P10	53.200-P10	



TRANSFER

17.056	17.017	17.018
Transfers include two screws (for open and closed tray).		




DIGITAL ANALOGUE

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

NP	SP	WP	Multi-Unit
3.3	3.75-4.2	5.0-6.0	4.8
14.056/D	14.095/D	14.096/D	14.007/D ¹
			14.085/D ²
			¹ Standard Multi-unit ² MIS® Multi-unit



ANALOGUE

10 Pack.


¹ Standard Multi-Unit ² MIS® Multi-Unit

14.056	14.095	14.096	14.007 ¹
			14.085 ²
	14.095-P10	14.096-P10	14.007-P10
			¹ Standard Multi-Unit ² MIS® Multi-Unit



TEMPORARY ABUTMENT

	24.095	24.096
	25.095	25.096



DESSLoc®

GH 1,0mm	88.017/1	88.018/1
GH 2,5mm	88.017/2	88.018/2
GH 3,5mm	88.017/3	88.018/3
GH 4,5mm	88.017/4	88.018/4
GH 5,5mm	88.017/5	




			NP	SP	WP	Multi-Unit
			3.3	3.75-4.2	5.0-6.0	4.8
TI-BASE						
	GH 0,3mm				15.096	15.085 <i>mis</i>
	GH 0,8mm		15.056	15.095		
	GH 1,5mm		15.256	15.295	15.296	
	GH 3,0mm		15.356	15.395	15.396	
	GH 0,3mm				16.096	
	GH 0,8mm		16.056	16.095		
	GH 1,5mm		16.256	16.295	16.296	
	GH 3,0mm		16.356	16.395	16.396	
TI-BASE CASTABLES						
	5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5
	5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5

MULTI-UNIT
ABUTMENTS

GH 1,0mm	0°			40.019/1	40.020/1
GH 2,0mm	0°			40.019/2	40.020/2
GH 3,0mm	0°			40.019/3	40.020/3
GH 4,0mm	0°			40.019/4	
GH 5,0mm	0°			40.019/5	
GH 2,5mm	17°			41.719/2	41.720/2
GH 3,5mm	17°			41.719/3	41.720/3
GH 3,5mm	30°			43.019/3	43.020/3
GH 4,5mm	30°			43.019/4	43.020/4

4.8
MULTI-UNIT PLATFORM



			NP	SP	WP	Multi-Unit
			3.3	3.75-4.2	5.0-6.0	4.8
AURUMBase®						
				36.017	36.018	
				37.017	37.018	
PLASTIC SLEEVES FOR AURUMBase®						
	5 Pack.	0°		33.100-P5	33.100-P5	
	5 Pack.	10°		33.101-P5	33.101-P5	
	5 Pack.	20°		33.102-P5	33.102-P5	



Torx® Screw AURUMBase®				
(replacement screw)			19.417	19.417
Recommended torque			30 Ncm	30 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



ABUTMENTS				
	0°		13.095	13.096
	15°		22.095	22.096
	25°		23.095	23.096



		NP	SP	WP	Multi-Unit
SCREW HEX 1,27mm		3.3	3.75-4.2	5.0-6.0	4.8
TIN		19.056	19.081	19.081	19.110 
		19.656	19.681	19.681	
	Recommended torque	30 Ncm	30 Ncm	30 Ncm	30 Ncm
TIN					19.085 
					19.685 
	Recommended torque				15 Ncm



PRE-MILLED BLANK					
Ø 10mm	Ti	61.056	61.095	61.096	
Ø 14mm	CoCr		71.095	71.096	




CoCr BASE					
		30.056	30.095	30.096	
		31.056	31.095	31.096	




CASTABLE					
10 Pack.		28.056	28.095	28.096	
			11.017-P10	11.018-P10	
		29.056	29.095	29.096	
10 Pack.			12.017-P10	12.018-P10	

AURUMBase®

Our angulated screw channel solution, with an angular adjustment of up to 25° and 360° freedom, will allow you to solve the most complex cases. A gold anodized finish with SelectGrip® enhanced surface treatment and a 3mm chimney will create the best aesthetic results without sacrifices.





SCAN
ABUTMENTS



		NP	SP	WP	Multi-Unit
					4.8
		52.081	52.082	52.083	52.085 <i>MIS</i>
10 Pack.		53.100-P10	53.100-P10	53.100-P10	
10 Pack.		53.200-P10	53.200-P10	53.200-P10	



TRANSFER

17.081 17.082 17.083

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*

14.081/D 14.082/D 14.083/D 14.007/D ¹
14.085/D ²

**Other retention screws available, see page 59*

¹ Standard Multi-Unit

² MIS® Multi-Unit



Platform Colour Coded

14.081 14.082 14.083 14.007 ¹
14.085 ²

¹ Standard Multi-Unit

² MIS® Multi-Unit



TI-BASE

		NP	SP	WP	Multi-Unit
					4.8
GH 0,3mm					15.085 <i>MIS</i>
GH 1,0mm		15.081	15.082	15.083	
GH 0,3mm					16.085 <i>MIS</i>
GH 1,0mm		16.081	16.082	16.083	



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.007-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	



SCREW HEX 1,27mm

	19.056	19.081	19.081	
				19.085 <i>MIS</i>
TIN	19.656	19.681	19.681	
TIN				19.685 <i>MIS</i>


Recommended torque

30 Ncm 30 Ncm 30 Ncm 25 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.081	61.082	61.083
--------	----	--------	--------	--------



CONIC MV

compatible with
Mis® V3

-  non-engaging
-  engaging



SCAN ABUTMENTS



- 10 Pack. 
- 10 Pack. 

	NP	SP	Multi-Unit
			4.8
	52.084	52.082	52.085 <i>MIS</i>
	53.100-P10	53.100-P10	
	53.200-P10	53.200-P10	




TRANSFER

17.084 17.082

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE



Retention screw Hex 1,27 mm*

14.084/D	14.082/D	14.007/D ¹
		14.085/D ²


**Other retention screws available, see page 59*

¹ Standard Multi-Unit

² MIS® Multi-Unit



ANALOGUE



14.084	14.082	14.007 ¹
		14.085 ²

¹ Standard Multi-Unit

² MIS® Multi-Unit



TI-BASE

		NP	SP	Multi-Unit
				4.8
GH 0,3mm				15.085 <i>MIS</i>
GH 1,0mm		15.084	15.082	
GH 1,0mm		16.084	16.082	



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.007-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	



SCREW HEX 1,27mm

	19.056	19.081	
			19.085
TIN	19.656	19.681	
TIN			19.685


Recommended torque

30 Ncm 30 Ncm 25 Ncm



PRE-MILLED BLANK


Ø 10mm	Ti	61.084	61.082
--------	----	--------	--------



CONIC ANYR

compatible with
Megagen AnyRidge®

-  non-engaging
-  engaging



SCAN
ABUTMENTS



50.057



52.057



10 Pack. 

53.100-P10



10 Pack. 

53.200-P10



TRANSFER

17.057

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*

14.057/D

**Other retention screws available, see page 59*



ANALOGUE



Platform Colour Coded

14.057



TI-BASE

GH 2,0mm



15.057

GH 3,0mm



15.357

GH 2,0mm



16.057

GH 3,0mm



16.357



TI-BASE CASTABLES

5 Pack.




32.001-P5

5 Pack.



33.001-P5



AURUMBase®



36.057



37.057

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.

0°

33.100-P5

5 Pack.

10°

33.101-P5

5 Pack.

20°

33.102-P5



Torx® Screw AURUMBase®

(replacement screw)

19.460

Recommended torque

35 Ncm

Torx® BALL SCREWDRIVER

20mm

24mm

30mm

35mm



DT20TB06

DT24TB06

DT30TB06

DT35TB06



SCREW HEX 1,20mm

		19.057
TIN		19.657
Recommended torque		35 Ncm




PRE-MILLED BLANK


Ø 10mm	Ti	61.057
Ø 14mm	CoCr	71.057



Tight spaces not a problem

Restorations where the interdental space is limited have always been challenging. Until now, existing solutions would have meant making sacrifices like reducing the amount of restorative material or achieving a poor aesthetic result. With a narrow shape and a 20 degree entry angle, ELLIPTIBase® is the most complete solution for these type of situations.





SCAN
ABUTMENTS



10 Pack. ○
10 Pack. △

52.077
53.100-P10
53.200-P10



TRANSFER

17.077

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

Platform Colour Coded

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

14.077/D



ANALOGUE

Platform Colour Coded

14.077



TI-BASE

GH 2,9mm ○ 15.077
GH 2,9mm ○ 16.077



TI-BASE CASTABLES

5 Pack. ○ 32.001-P5
5 Pack. △ 33.001-P5



SCREW HEX 1,18mm

19.077
TIN 19.677

Recommended torque 25 Ncm

On Abutment® 19.062


TIN 19.662

Recommended torque 15 Ncm



PRE-MILLED BLANK

Ø 10mm Ti 61.077



SCAN
ABUTMENTS



10 Pack. 

10 Pack. 

50.080

52.080

53.100-P10

53.200-P10



TRANSFER

17.080

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

 Platform Colour Coded

Retention screw Hex 1,27 mm*

*Other retention screws available, see page 59

14.080/D



ANALOGUE

 Platform Colour Coded

14.080



TI-BASE



15.080



16.080



TI-BASE CASTABLES

5 Pack.




32.001-P5

5 Pack.



33.001-P5



AURUMBase®

GH 2,2mm



36.080

GH 0,9mm



36.180

GH 2,2mm



37.080

GH 0,9mm



37.180

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.

0°

33.100-P5

5 Pack.

10°

33.101-P5

5 Pack.

20°

33.102-P5



Torx® Screw AURUMBase®

2,2mm

19.444

(replacement screw)

0,9mm

19.443

Recommended torque

30 Ncm

Torx® BALL SCREWDRIVER

20mm

24mm

30mm

35mm

DT20TB06

DT24TB06

DT30TB06

DT35TB06





SCREW HEX 1,40mm

		19.080
	TIN	19.680
Recommended torque		30 Ncm




PRE-MILLED BLANK

Ø 10mm	Ti	61.080
Ø 14mm	Ti	62.080
Ø 14mm	CoCr	71.080

The importance of a precise image

Our Ball desktop scan abutment designed in conjunction with Exocad® and the revolutionary toolless intraoral scan abutment, have been carefully designed to produce the most precise image. Intraoral scan abutment manufactured in titanium allows for Rx control, now with ZrN coating to reduce glitter effect.





SCAN
ABUTMENTS



10 Pack. 

10 Pack. 

52.086

53.100-P10

53.200-P10



TRANSFER

17.086

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*

14.086/D

**Other retention screws available, see page 59*



ANALOGUE



Platform Colour Coded

14.086



TI-BASE



15.086



16.086



TI-BASE CASTABLES

5 Pack.



32.001-P5

5 Pack.



33.001-P5



SCREW HEX 1,22mm

19.086

TIN

19.686

On Abutment

19.035

19.036

TIN


19.635

TIN

19.636

Recommended torque


15 Ncm



CONIC DENT

compatible with
Dentium Superline™ & Implantium®

- non-engaging
- ⬡ engaging



SCAN
ABUTMENTS



- 10 Pack. ○
- 10 Pack. ⬡

RP				
3.6	3.8	4.0	4.5	5.0
50.090				
52.090				
53.100-P10				
53.200-P10				



TRANSFER

17.090

Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

Platform Colour Coded

Retention screw Hex 1,27 mm*

14.090/D

*Other retention screws available, see page 59



Platform Colour Coded

ANALOGUE

RP				
3.6	3.8	4.0	4.5	5.0
14.090				



TI-BASE

GH 0,9mm	○	15.090
GH 1,5mm	○	15.190
GH 3,0mm	○	15.390
GH 0,9mm	⬡	16.090
GH 1,5mm	⬡	16.090
GH 3,0mm	⬡	16.390




TI-BASE CASTABLES

5 Pack.	○	32.001-P5
5 Pack.	⬡	33.001-P5



			RP				
			3.6	3.8	4.0	4.5	5.0
MULTI-UNIT ABUTMENTS							
GH 1,5mm	0°	◯	40.090/1				
GH 2,5mm	0°	◯	40.090/2				
GH 3,5mm	0°	◯	40.090/3				
GH 2,5mm	17°	◯	41.790/2				
GH 3,5mm	17°	◯	41.790/3				
GH 2,5mm	30°	◯	43.090/2				
GH 3,5mm	17°	◯	43.090/3				



AURUMBase®		
	◯	36.090
	◯	37.090



PLASTIC SLEEVES FOR AURUMBase®		
5 Pack.	0°	33.100-P5
5 Pack.	10°	33.101-P5
5 Pack.	20°	33.102-P5



Torx® Screw AURUMBase®	
(replacement screw)	19.405
Recommended torque	30 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06



		RP				
		3.6	3.8	4.0	4.5	5.0
SCREW HEX 1,28mm						
	TITAN	19.090				
	TITANium	19.690				
Recommended torque		30 Ncm				
On Abutment®		19.066				
	TITAN	19.666				
Recommended torque		10 Ncm				




PRE-MILLED BLANK		
Ø 10mm	Ti	61.090
Ø 14mm	Ti	62.090
Ø 14mm	CoCr	71.090



CoCr BASE	
◯	30.090
◯	31.090




CASTABLE	
◯	28.090
◯	29.090



INTERNAL SEM

compatible with
Bego Semados®

- non-engaging
- engaging




SCAN ABUTMENTS

			3.25/3.75	4.1	4.5	5.5
			52.091	52.092	52.093	52.094
	10 Pack.		53.100-P10	53.100-P10	53.100-P10	53.100-P10
	10 Pack.		53.200-P10	53.200-P10	53.200-P10	53.200-P10



TRANSFER

			17.091	17.092	17.093	17.094
Transfers include two screws (for open and closed tray).						



DIGITAL ANALOGUE

	Platform Colour Coded	Retention screw Hex 1,27 mm*	14.091/D	14.092/D	14.093/D	14.094/D
*Other retention screws available, see page 59						



ANALOGUE

	Platform Colour Coded		14.091	14.092	14.093	14.094
--	-----------------------	--	--------	--------	--------	--------



			3.25/3.75	4.1	4.5	5.5
GH 0,7mm					15.091	15.094
GH 0,9mm				15.091		
GH 1,1mm		15.091				
GH 0,7mm					16.091	16.094
GH 0,9mm				16.091		
GH 1,1mm		16.091				



TI-BASE CASTABLES

5 Pack.		32.001-P5	32.001-P5	32.001-P5	32.001-P5
5 Pack.		33.001-P5	33.001-P5	33.001-P5	33.001-P5




SCREW HEX 1,27mm

		19.091	19.091	19.091	19.091
TIN		19.691	19.691	19.691	19.691
Recommended torque		30 Ncm	30 Ncm	30 Ncm	30 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.091	61.092	61.093	61.094
Ø 14mm	Ti				62.094
Ø 14mm	CoCr				71.094



CONIC BIO
compatible with
Biotech KONTACT®

- non-engaging
- ⚙️ engaging



SCAN
ABUTMENTS

		NP	RP
		52.078	
	10 Pack. ○	53.100-P10	
	10 Pack. ⚙️	53.200-P10	



TRANSFER

17.078
Transfers include two screws (for open and closed tray).



DIGITAL
ANALOGUE

	Platform Colour Coded	Retention screw Hex 1,27 mm*	14.078/D
*Other retention screws available, see page 59			



ANALOGUE

	Platform Colour Coded	14.078
--	-----------------------	--------



TI-BASE

		NP	RP
	○	15.078	
	⚙️	16.078	



TI-BASE CASTABLES

5 Pack.	○	32.001-P5
5 Pack.	⚙️	33.001-P5



SCREW HEX 1,27mm

	19.078	19.063	19.064
TIN	19.678	19.663	19.664
Recommended torque	20 Ncm	20 Ncm	20 Ncm



PRE-MILLED BLANK

Ø 10mm	Ti	61.078
Ø 14mm	CoCr	71.078



INTERNAL BT

compatible with
BTI®

- non-engaging
- engaging



SCAN ABUTMENTS



- 10 Pack.
- 10 Pack.

4.1	MULTI IM®
52.079	
53.100-P10	
53.200-P10	



TRANSFER

17.079

Transfers include two screws (for open and closed tray).



DIGITAL ANALOGUE

Platform Colour Coded

Retention screw Hex 1,27 mm*
*Other retention screws available, see page 59


14.079/D



ANALOGUE

Platform Colour Coded

14.079



TI-BASE

-
-

4.1	MULTI IM®
15.079	
16.079	



TI-BASE CASTABLES

- 5 Pack.
- 5 Pack.

32.001-P5
33.001-P5



SCREW HEX 1,27mm

TIN

19.079	19.065
19.679	19.665

Recommended torque


35 Ncm	10 Ncm
--------	--------

EXTERNAL
HEX PH

compatible with
PHIBO® TSH®

 non-engaging

 engaging



SCAN
ABUTMENTS

		NP	SP	WP
		S2	S3-S4	S5
		52.030	52.031	52.032
	10 Pack. 	53.100-P10	53.100-P10	53.100-P10
	10 Pack. 	53.200-P10	53.200-P10	53.200-P10



TRANSFER

		NP	SP	WP
		17.030	17.031	17.032
Transfers include two screws (for open and closed tray).				



DIGITAL
ANALOGUE




		14.030/D	14.031/D	14.032/D
Retention screw Hex 1,27 mm* *Other retention screws available, see page 59				



ANALOGUE



		14.030	14.031	14.032
--	--	--------	--------	--------



		NP	SP	WP
		S2	S3-S4	S5
		36.030	36.031	36.032
		36.030/2	36.031/2	36.032/2
		36.030/3	36.031/3	36.032/3
		37.030	37.031	37.032
		37.030/2	37.031/2	37.032/2
		37.030/3	37.031/3	37.032/3

PLASTIC SLEEVES FOR AURUMBase®

	5 Pack.	0°	33.100-P5	33.100-P5	33.100-P5
	5 Pack.	10°	33.101-P5	33.101-P5	33.101-P5
	5 Pack.	20°	33.102-P5	33.102-P5	33.102-P5



Torx® Screw AURUMBase®

(replacement screw)	19.403	19.402	19.402
Recommended torque	35 Ncm	35 Ncm	35 Ncm



Torx® BALL SCREWDRIVER	20mm	24mm	30mm	35mm
	DT20TB06	DT24TB06	DT30TB06	DT35TB06


MULTI-UNIT
CONNECTION

 non-engaging



HEALING
ABUTMENT

			RP MULTI-UNIT 4.8	WP MULTI-UNIT 6.5
		Uni Grip	20.100	20.054
		Hex 1,20mm	20.100H120	
		Hex 1,22mm	20.100H122	
		Hex 1,27mm	20.100H127	
		Torx® Ball	20.100TB	



SCAN
ABUTMENTS

			50.007	50.047
			52.007	52.047
	5 Pack.		53.107-P5	53.154-P5
	5 Pack.		53.107-P5	53.155-P5



TRANSFER

			17.007	17.050
5,0mm			17.507	
5,0mm	4 Pack.		17.507-P4	
8,0mm			17.807	
8,0mm	4 Pack.		17.807-P4	



DIGITAL
ANALOGUE



Platform Colour Coded

Retention screw Hex 1,27 mm*

**Other retention screws available, see page 59*



ANALOGUE



Platform Colour Coded

10 Pack.



TEMPORARY
ABUTMENT

	 Uni Grip	24.007	24.030
	 Hex 1,20mm	24.007H120 19.055	
	 Hex 1,22mm	24.007H122 19.065	
	 Hex 1,27mm	24.007H127 19.110	
	 Torx® Ball	24.007TB 19.306	



TEMPORARY
ABUTMENT PEEK

	 Uni Grip	25.707	25.754
	 Hex 1,20mm	25.707H120 19.055	
	 Hex 1,22mm	25.707H122 19.065	
	 Hex 1,27mm	25.707H127 19.110	
	 Torx® Ball	25.707TB 19.306	



DESSLoc®

RP
MULTI-UNIT
4.8

WP
MULTI-UNIT
6.5

88.007

TI-BASE



6,0mm		Uni Grip	15.007	15.054
6,0mm		Hex 1,20mm	15.007H120 19.055	
6,0mm		Hex 1,22mm	15.007H122 19.065	
6,0mm		Hex 1,27mm	15.007H127 19.110	
6,0mm		Torx® Ball	15.007TB 19.306	
8,0mm			15.007/8	

TI-BASE CASTABLES



5 Pack.		32.007-P5	32.054-P5
5 Pack.		32.007/8-P5	

AURUMBase®



	36.007	36.054
---	--------	--------

PLASTIC SLEEVES FOR AURUMBase®



5 Pack.	0°	33.700-P5	33.700-P5
5 Pack.	10°	33.701-P5	33.701-P5
5 Pack.	20°	33.702-P5	33.702-P5

Torx® Screw AURUMBase®



(replacement screw)	19.306	19.331
Recommended torque	15 Ncm	15 Ncm

Torx® BALL SCREWDRIVER



20mm	24mm	30mm	35mm
DT20TB06	DT24TB06	DT30TB06	DT35TB06



SCREW UG

RP
MULTI-UNIT
4.8

WP
MULTI-UNIT
6.5

	19.006	19.031
10 Pack.	19.006-P10	19.031-P10
DLC	19.206	19.231
TIN	19.606	19.631



	Hex 1,20mm	19.055
	Hex 1,22mm	19.065
	Hex 1,27mm	19.110

Recommended torque

15 Ncm 15 Ncm

CoCr BASE



6,0mm		30.007	30.054
6,0mm		Hex 1,20mm	30.007H120 19.055
6,0mm		Hex 1,22mm	30.007H122 19.065
6,0mm		Hex 1,27mm	30.007H127 19.110
6,0mm		Torx® Ball	30.007TB 19.306
For wealding techniques		30.107	

CASTABLE



		28.007	28.054
10 Pack.		11.007-P10	11.054-P10



SCAN
ABUTMENTS

10 Pack.

53.153-P10



TRANSFER

17.054




DIGITAL
ANALOGUE

Retention screw Hex 1,27 mm*

14.054/D

**Other retention screws available, see page 59*

ANALOGUE



10 Pack.

14.054

14.054-P10



0° - 20° x1
0° - 20° x2


DLRK-P1
DLRK-P2



0° - 40° x1
0° - 40° x2

DLER-P1
DLER-P2

RETENTION KITS



MALES

0° - 20° x4
0° - 20° x4
0° - 20° x4
0° - 40° x4
0° - 40° x4
0° - 40° x4

DLN2-P4
DLN3-P4
DLN4-P4
DLN1-P4
DLN6-P4
DLN7-P4



SPECIAL TOOLS

Direct to torque wrench Nobel Biocare™

DTACTLOC



Direct to torque wrench Straumann®

DTSTRLOC



DESSLoc®

DTLOC



Adaptor from Straumann® torque wrench to square connection

DTSTSQ




Nobel Biocare™ Adapter from 8mm to square connection




DTNBSQ

DESSLoc[®] retention elements removal and insertion tool

EXTRACTION



1. Place the tool in the following position, making two complete turns on the gold side.



2. Use the tool to remove the old retainer. Press with it until you hear a “click”, finish with a slight lever movement.
3. Screw in the tool to eject the retainer.

RETENTION INSERTION



4.



5.




6.

4. Completely remove the gold part of the tool.
5. Place the new retainer on the end.
6. Press on the cup until you hear a “click”.

DESSLoc[®] INSERTION



1. Use the back of the PHLOC tool to facilitate placement on the models. Fit the tool with the upper part of the DESSLoc[®]



PURE SWITCH^{by} DESS[®]

“simple switch, perfect connection”

The Pure Switch concept is the philosophy upon which our products are conceived. The best example of this are our screws which have been manufactured using the same dimensions and characteristics as the implant branded screws. As the key component of every ensemble, DESS[®] screws have a perfect fit and can be interchanged with the branded version without any compromise.



HOLDERS



HLDSTD
For standard
98mm disc machines

Any milling machine with
a 98mm disc holder



HLDIMES
For Imes-Icore® machines

Coritec 150i pro
Coritec 245i touch
Coritec 250i touch
Coritec 350i / Coritec 250i Loader
Coritec 450i
Coritec 550i
Coritec 650i / Coritec 650i Loader



HLDVHFR5
For VHF®

VHF® R5




HLDYENA
For Yenadent®

D6
D15
D43N



HLDVHFZ4
For VHF®

VHF® Z4



HLD3R
For machines with
support / changer 3R

D14P
DC40
D43N



HLDVHF
For VHF®

VHF® S1
VHF® S2
VHF® K5
VHF® S5




HLDVHFN4
For VHF®

VHF® N4

AURUM MILLING TOOL



	ØA	ØB	ØC	L1	L2	L3	R
MT3D16L5	3	1,6	0,95	5	51	0,5	0,1
MT4D16L5	4	1,6	0,95	5	51	0,5	0,1
MT6D16L5	6	1,6	0,95	5	51	0,5	0,1



BLOCKS




Ti	□ 15mm	L 20mm	63.020
Ti	□ 15mm	L 40mm	63.040
CoCr	□ 15mm	L 20mm	73.020
CoCr	□ 15mm	L 40mm	73.040


TOOLS

SCREWDRIVERS TIPS

HEX	18mm	20mm	24mm	30mm	35mm	compatible with
0,90mm		DTC20H090	DTC24H090			Neodent® CM
1,00mm		DTC20H100	DTC24H100		DTC35H100	Dentsply Ankylos® C/X
1,20mm	DTC18H120	DTC20H120	DTC24H120	DTC30H120	DTC35H120	3i Osseotite®, Neodent® CM & GM, 3iCertain®, Megagen AnyRidge®
1,22mm		DTC20H122	DTC24H122		DTC35H122	Dentsply Friadent Xive®, Osstem® TS, GlobalD®, BTI®, Biotech Kontakt®
1,25mm		DTC20H125	DTC24H125			
1,27mm		DTC20H127	DTC24H127	DTC30H127	DTC35H127	Astra tech™, OsseosSpeed™, Zimmer Screw-Vent® & Mis® Seven/ C1 & V3, BioHorizons®, AstraTech Implants System EV™, Camlog®, Bego Semados®
1,30mm		DTC20H130	DTC24H130			
1,50mm		DTC20H150	DTC24H150			
1,70mm		DTC20H170	DTC24H170			
1,75mm		DTC20H175	DTC24H175			



Torx® 06






Torx® ball



DT20T06	DT24T06	DT30T06	DT35T06	Straumann® Soft Tissue Level & synOcta® Straumann® Bone Level
DT20TC06	DT24TC06	DT30TC06	DT35TC06	

DT20TB06	DT24TB06	DT30TB06	DT35TB06	AURUMBase®, ELLIPTIBase®
----------	----------	----------	----------	--------------------------



UG	18mm	20mm	24mm	30mm	35mm	compatible with
		DT20U	DT24U	DT30U	DT35U	Nobel Brånemark®, NobelReplace® Select™, NobelActive™ & NobelReplace® CC, Multi-unit®
NEO						Grand Morse® Neodent®
			DT24NG	DT30NG	DT35NG	
Squared						
1,30mm		DT20S130	DT24S130			
						

SCREWDRIVER HANDLES

Premium laboratory handle



PHLABGOLD



PHLABBLACK

Dentist handle



7mm

MDTMAN700



9mm

MDTMAN900

DESSLoc® Removal & Insertion Tool



PHLOC

TI-BASE EXTRACTING TOOL

	For Conelog® NP/RP	DCCON
	For Conelog® WP	DCCON50
	For ZIMMER SCREW-VENT®	DCZIM
	For Medentis ICX®	DCICX

TORQUE WRENCH TOOL

	Type Straumann 7mm	DTSTTW
	Type Nobel 8mm	DTNBTW

TRAY MULTI-USE

	BOX1
--	------

SPECIAL TOOLS

	Multi-unit® RP	DTMURP
	Multi-unit® WP	DTMUWP
	UniAbutment®	DTUNI
	Octagon	DTSYN
	Torque wrench adaptor	TWADAP
	DESSLoc®	DTLOC
	Adaptor from 7mm to Straumann® torque wrench to square connection	DTSTSQ
	Nobel Biocare™ adapter from 8mm to square connection	DTSTWH
	Direct to torque wrench Nobel Biocare™ 7mm	DTACTLOC
	Direct to torque wrench Straumann® 8mm	DTSTRLOC
	Straumann® torque wrench	DTSTCA
	Screwdriver tips adaptor from Straumann® torque wrench	DTNBSQ
	Screwdriver tips adaptor from Nobel Biocare™ torque wrench	DTNBWH
	Nobel Biocare™ contra-angles adapter	DTNBCA
	Handpiece adaptor for Square connection	DTANGSQ
	Manual driver for Intraoral Scan abutment	DTPEI
	Manual tool for Intraoral Scan abutment	DTPEIPEEK-P2



CLINICAL SCREWS RECOMMENDED TORQUE

NOBELACTIVE® &
NOBELREPLACE® CC
ACTIVE HEX



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
3.0	19.068	15	Nobel Biocare® 36805
3.0	19.268	15	Nobel Biocare® 37890
3.0	19.668	15	
NP	19.027	35	Nobel Biocare® 31168
NP	19.227	35	Nobel Biocare® 29282
NP	19.101	35	
NP	19.327	35	
NP	19.627	35	
NP AURUMBase®	19.441	35	
RP	19.030	35	Nobel Biocare® 28815
RP	19.230	35	Nobel Biocare® 36918
RP	19.328	35	
RP	19.630	35	
RP AURUMBase®	19.442	35	
On1®	19.021	35	Nobel Biocare® 38751
On1®	19.221	35	Nobel Biocare® 38750
On1®	19.621	35	
NP on implant for Multi-unit®	19.507	35	
RP on implant for Multi-unit®	19.508	35	

NP	19.004	35	Nobel Biocare® 31170
NP	19.204	35	Nobel Biocare® 36818
NP	19.304	35	
NP	19.604	35	
NP AURUMBase®	19.404	35	
RP/WP/6.0	19.005	35	Nobel Biocare® 29293
RP/WP/6.0	19.205	35	Nobel Biocare® 29475
RP/WP/6.0	19.305	35	
RP/WP/6.0	19.605	35	
RP/WP AURUMBase®	19.405	35	
NP on implant for Multi-unit®	19.504	35	
RP on implant for Multi-unit®	19.505	35	

NOBELREPLACE
SELECT™
TRI-LOBE



NOBEL BRÅNEMARK®
EXT HEX UNIVERSAL



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
NP	19.027	35	Nobel Biocare® 31168
NP	19.227	35	Nobel Biocare® 29282
NP	19.327	35	
NP	19.101	35	
NP	19.627	35	
NP AURUMBase®	19.401	35	
RP	19.028	35	Nobel Biocare® 29290
RP	19.228	35	Nobel Biocare® 29283
RP	19.328	35	
RP	19.628	35	
RP AURUMBase®	19.402	35	
WP	19.029	35	Nobel Biocare® 31169
WP	19.229	35	Nobel Biocare® 29284
WP	19.329	35	
WP	19.629	35	
NP on implant for Multi-unit®	19.501	35	
RP on implant for Multi-unit®	19.502	35	

NP/RP/WP	19.034	20	Biomet 3i® UNIHT
NP/RP/WP	19.634	20	
NP/RP/WP AURUMBase®	19.402	20	
RP on implant for Multi-unit®	19.512	20	
Low Profile® Abutment ø4.8	19.106	10	Biomet 3i® LPCTSH

NP/RP/WP	19.011	20	Biomet 3i® IUNIHT
NP/RP/WP	19.311	20	
NP/RP/WP	19.611	20	
NP/RP/WP AURUMBase®	19.438	20	
Low Profile® Abutment ø4.8	19.106	10	Biomet 3i® LPCTSH

NNC	19.070	35	Straumann® AG 048.314
NNC	19.670	35	
RN/WN on implant	19.026	35	Straumann® AG 048.356
RN/WN on implant	19.102	35	Straumann® AG 049.181
RN/WN on implant	19.328	35	
RN/WN on implant	19.626	35	
RN/WN AURUMBase®	19.446	35	
RN/WN on synOcta® abutment	19.009	15	Straumann® AG 048.350
RN/WN on synOcta® abutment	19.103	15	
RN/WN on synOcta® abutment	19.609	15	

3i CERTAIN®
EXT HEX USA



3i OSSEOTITE®
INT HEX CLICK



STRAUMANN® SOFT
TISSUE LEVEL &
SynOcta®
OCTAGON



STRAUMANN®
BONE LEVEL
CONICAL BL



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
NC on implant	19.032	35	Straumann® AG 025.2908
NC on implant	19.632	35	
RC on implant	19.033	35	
RC on implant	19.633	35	
NC/RC on implant	19.037	35	Straumann® AG 025.4900
NC/RC on implant	19.104	35	
NC/RC on implant	19.311	20	Straumann® AG 025.4906
NC/RC AURUMBase®	19.443	35	
NC/RC on implant for Multi-unit®	19.544	35	
NC/RC on Multi-base® abutment	19.041	15	
NC/RC on Multi-base® abutment	19.105	15	Straumann® AG 025.0900-04
NC/RC on SRA® abutment	19.098	15	

STRAUMANN®
BLX
CONICAL BLX



RB/WB	19.054	35	
RB/WB	19.700	35	
RB/WB AURUMBase®	19.441	35	

NEODENT®
GRAND MORSE®
NEO GM



GM	19.051	20	
GM AURUMBase®	19.444	20	

ZIMMER
SCREW-VENT®
INT HEX USA



NP/RP/WP	19.012	30	Zimmer® URS2
NP/RP/WP	19.012L	30	
NP/RP/WP	19.042	30	Zimmer® MHLAS
NP/RP/WP	19.642	30	
NP/RP/WP	19.312	30	
NP/RP/WP AURUMBase®	19.417	30	
NP (MIS®)	19.327	35	
NP on implant for Multi-unit®	19.519	30	
RP on Multi-unit® abutment ø4.8	19.110	15	
RP on Multi-unit® abutment ø4.8	19.300	15	

ASTRA TECH
IMPLANT
SYSTEM™ EV
CONIC EVO



PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
3.6 on implant	19.059	25	
3.6 on implant	19.659	25	
3.6 AURUMBase®	19.441	25	
4.2 on implant	19.060	25	
4.2 on implant	19.660	25	Astra Tech® 25.205
4.2 AURUMBase®	19.460	25	
4.8 on implant	19.061	25	
4.8 on implant	19.661	25	
4.8 AURUMBase®	19.461	25	
UniAbutment® EV	19.058	15	
UniAbutment® EV	19.658	15	

ASTRA TECH
OSSEOSPEED™
INT HEX CONIC



NP on implant	19.015	15	Dentsply® 24189
NP on implant	19.615	15	
RP on implant	19.016	20	Dentsply® 24449
RP on implant	19.616	20	
RP AURUMBase®	19.441	20	
WP on implant	19.017	25	Dentsply® 24209
WP on implant	19.617	25	
WP AURUMBase®	19.425	25	
UniAbutment® 20°/45°	19.018	15	Dentsply® 22435
UniAbutment® 20°/45°	19.618	15	

DENTSPLY
FRIADENT XIVE®
INT HEX FD



NP/RP/WP	19.024	24	Dentsply® 454305
NP/RP/WP	19.624	24	
NP/RP/WP AURUMBase®	19.438	24	

DENTSPLY
ANKYLOS® C/X
INT ANK



All platforms	19.048	15	
All platforms	19.050	15	
All platforms	19.648	15	

	PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
OSSTEM® TS CONIC OSS 	Mini on implant	19.072	20	Osstem® GSABSST
	Mini on implant	19.672	20	
	Mini AURUMBase®	19.444	20	
	Standard on implant	19.073	30	
	Standard on implant	19.673	30	
	Standard AURUMBase®	19.461	30	
CAMLOG® INT CAM 	3,3 / 3,8 / 4,3 on implant	19.044	20	Camlog® J4005. 1601
	3,3 / 3,8 / 4,3 on implant	19.644	20	
	3,8 / 4,3 AURUMBase®	19.444	20	Camlog® J4005. 2001
	5.0 on implant	19.045	20	
	5.0 on implant	19.645	20	
CONELOG® INT CON 	3,3 / 3,8 / 4,3 on implant	19.074	20	Camlog® C4005. 1601
	3,3 / 3,8 / 4,3 on implant	19.674	20	Camlog® C4005. 2001
	5.0 on implant	19.076	20	
	5.0 on implant	19.676	20	
BIOHORIZONS® INTERNAL INT HEX BH 	On implant	19.043	30	BioHorizons® PXAS
	On implant	19.043	30	Nobel Biocare® 29287
	Multi-unit® RP ø4.8	19.006	15	
	Low Profile® abutment ø4.8	19.106	10	Biomet 3i® LPCTSH
	Multi-unit® RP ø4.8	19.110	15	Nobel Biocare® 29285
	Multi-unit® RP ø4.8	19.206	15	
	Multi-unit® RP ø4.8 AURUMBase®	19.306	15	
	Multi-unit® RP ø4.8	19.606	15	Nobel Biocare® 31163
	MIS® Multi-unit® RP ø4.8	19.300	15	
	Multi-unit® WP ø6.0	19.031	15	
	Multi-unit® WP ø6.0	19.231	15	Nobel Biocare® 29286
	Multi-unit® WP ø6.0 AURUMBase®	19.331	15	
	Multi-unit® WP ø6.0	19.631	15	

	PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
BIOHORIZONS® EXTERNAL EXT HEX BH 	NP	19.013	30	BioHorizons® 130-050
	NP	19.312	30	BioHorizons® 140-050
	RP/WP	19.014	30	
	RP/WP	19.328	35	
MIS® SEVEN INT HEX MI 	NP	19.056	30	MIS® MN-S0160
	NP	19.656	30	
	SP/WP	19.081	30	MIS® MD-S0220
	SP/WP	19.681	30	MIS® MU-S0220
	SP/WP AURUMBase®	19.417	30	
	Multi-unit® MIS®	19.085	25	
	Multi-unit® MIS®	19.685	25	MIS® MU-S0220
	Multi-unit® MIS®	19.685	25	
MIS® C1 CONIC MC 	NP	19.056	30	MIS® MN-S0160
	NP	19.656	30	
	SP/WP	19.081	30	MIS® MD-S0220
	SP/WP	19.681	30	MIS® MU-S0220
	Multi-unit® MIS®	19.085	25	
	Multi-unit® MIS®	19.685	25	
MIS® V3 CONIC MV 	NP	19.056	30	MIS® MN-S0160
	NP	19.656	30	
	SP	19.081	30	MIS® MD-S0220
	SP	19.681	30	MIS® MU-S0220
	Multi-unit® MIS®	19.085	25	
	Multi-unit® MIS®	19.685	25	
MEGAGEN AnyRidge® CONIC ANYR 	All platforms	19.057	35	Megagen® AANMSF
	All platforms	19.657	35	
	All platforms AURUMBase®	19.460	35	

	PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
ANTHOGYR AXIOM® CONIC AXIO 	On implant	19.077	25	Anthogyr® OPTS161
	On implant	19.677	25	
	Multi-unit® abutment	19.062	15	
	Multi-unit® abutment	19.662	15	

MEDENTIS ICX® CONIC IC 	On implant	19.080	30	Medentis® C-007-000001-01
	On implant	19.680	30	
	On implant AURUMBase® H2.2	19.444	30	
	On implant AURUMBase® H0.9	19.443	30	

GLOBALD IN-KONE® CONIC INK 	On implant	19.086	15	GlobalD® DVPCI
	On implant	19.686	15	
	Conical abutment Straight	19.035	15	GlobalD® DVPIINK
	Conical abutment Straight	19.635	15	
	Conical abutment angled	19.036	15	
	Conical abutment angled	19.636	15	

DENTIUM SUPERLINE™ & IMPLANTIUM® CONIC DENT 	On implant	19.090	30	Dentium® ASC 2045
	On implant	19.690	30	
	On abutment	19.066	10	Dentium® SRS 18T
	On abutment	19.666	10	

BEGO SEMADOS® INTERNAL SEM 	On implant	19.091	30	Bego® 57.031
	On implant	19.691	30	

	PLATFORM SIZE	DESS REF. N°	DESS RECOMMENDED TORQUE (Ncm)	ORIGINAL EQUIVALENT REF. N°
BIOTECH KONTACT® CONIC BIO 	On implant	19.078	20	Biotech® KVP
	On implant	19.678	20	
	Conical abutment Narrow	19.063	20	
	Conical abutment Narrow	19.663	20	
	Conical abutment Standard	19.064	20	Biotech® KECGMV
	Conical abutment Standard	19.664	20	

BTI® INTERNAL BT 	On implant	19.079	35	
	On implant	19.679	35	
	Multi IM®	19.065	10	
	Multi IM®	19.665	10	

PHIBO® TSH® EXTERNAL HEX PH 	NP AURUMBase®	19.403	35	
	SP/WP AURUMBase®	19.402	35	

MULTI-UNIT 	NP	19.006	15	
	NP	19.206	15	
	NP	19.606	15	
	WP	19.031	15	
	WP	19.231	15	
	WP	19.631	15	
	RP AURUMBase®	19.306	15	
	WP AURUMBase®	19.631	15	



DESS

> DENTAL SMART SOLUTIONS



DESS HQ

6.200m² Barberà del Vallès, Barcelona Spain



All DESS products are manufactured according to ISO 9001 and ISO 13485 with CE marking for all products
Class IIb

All DESS® AURUMBase® screws and screwdrivers are based on the standard ISO 10664

3i® and 3i Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA

ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XiVE®, FRIADENT®, ASTRATECH implant system™ EV and UniAbutment®

Cerec® are registered trademarks of DENTSPLY SIRONA Implants

BioHorizons® is a registered trademark of BioHorizons, INC.

Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and Multi-unit® are registered trademarks of Nobel Biocare AB,

Neodent®, Grand Morse® , Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland

MIS® is a registered trademark of MIS Implant Technologies Ltd.

Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA

Dyna® is a registered trademark of Dyna Dental Engineering BV

Locator® is a registered trademark of ZEST IP Holdings

Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG

Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD

TORX® is a trademark of ACUMENT INTELLECTUAL PROPERTIES, LLC

Bego Semados® is a registered trademark of BEGO GmbH & Co. KG

Anthogyr® is a registered trademark of Anthogyr SAS

Medentis ICX® is a registered trademark of medentis medical GmbH

Osstem®, Hiossen® is a registered trademark of OSSTEM IMPLANT CO. LTD

Dentium® is a trademark of Dentium Co. Ltd.

GOBALD® is a trademark of Global D SASU.

Biotech Kontakt® is a trademark Biotech Dental SAS.

BTI® is a trademark BTI Biotechnology Institute.

PHIBO® is a trademark PHIBO Dental Solutions S.L.

DESS products can be used only by licensed dentists, physicians and labs.
Please refer ever to the instructions of use available on the website
www.dessdental.com
or supplied with the product.

Some products may not be available in your country. Some of the products are cleared for US and Canada.
Please contact your local distributor for more details.



DOC-EN-05-FRC



www.dessdental.com